INSTRUCTIONS FOR USE

SINGLE USE NITRILE GLOVES

1140

FR	EN	DE	ES	IT
PT	SV	NL	PL	AR




CH	REP
-----------	------------

Medicom SAS
49124 Saint-Barthélemy-d'Anjou
France

Finagex SA
Rue Pierre-Fatio 12
1204 Genève
Switzerland

medicom-eu.com

SYMBOLS & PICTOGRAMS

FR

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EN

1140_MEDICOM KOLMI COVERFEEL TOUCH - INSTRUCTIONS FOR USE FOR CATEGORY III PERSONAL PROTECTIVE GLOVES, LETHAL OR IRREVERSIBLE RISKS - CLASS I MEDICAL EXAMINATION GLOVES - SINGLE USE - NON STERILE. This glove is recommended for protection against microorganisms (bacteria, fungi and viruses), tested chemicals and chemotherapy drugs. This glove is also intended for medical use (examination glove) by healthcare professionals. This glove is suitable for food contact. The glove is available in the following sizes: 6/XS - 7/S - 8/M - 9/L - 10/XL - 11/XL. **INSTRUCTIONS:** Choose the appropriate size. Put on glove(s) to the maximum. To achieve the desired protection, this glove can be used with other suitable PPE (coveralls, footwear, APR, ...). Be careful not to use non-recommended products. Pay attention to the symbols of danger.

EXPLANATION OF SYMBOLS, PICTOGRAMS AND GLOVE PERFORMANCE: (01) EN ISO 21420:2020 - Please read these instructions before using the products. (02) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (03) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (04) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test.

(05) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (06) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (07) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (08) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test.

(09) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (10) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (11) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (12) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test.

(13) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (14) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (15) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (16) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test.

(17) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (18) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (19) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (20) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test.

(21) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (22) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (23) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (24) EN ISO 374-5: 2016 - Protection against virus, bacteria and fungi. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test.

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