Medicom SAS
Bld de la chanterie | 49124 Saint Barthélemy d'Anjou
[BP 10059 – 49181 Saint Barthélemy d'Anjou Cedex]
France
Tel +33 0 241 963 434 | sales.support.pro.eu@medicom.com | www.medicom-eu.com

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ENGLISH – EN

This device is marked CE in accordance with Regulation (EU) 2017/745 on medical devices. Manufactured and packaged (in dispenser box) in France at Medicom SAS.

Composition: The mask is made from polypropylene, elastane and steel. Latex-, phthalate- and bisphenol-free.

Mask performance: This mask has been tested in accordance with the applicable standard EN14683:2019+AC:2019 and meets the following requirements according to the type indicated on the packaging:

	Type II	Type IIR
Bacterial filtration efficiency	≥98%	≥98%
Respirability	< 40 Pa/cm ²	< 60 Pa/cm ²
Splash resistance	Not applicable	16 kPa
Microbial cleanliness	≤30 cfu/g	≤30 cfu/g
Biocompatibility	non-cytotoxic (ISO 10993-5) non-irritant and non-sensitising (ISO 10993-10)	

Definition of label traceability pictograms :

Description of the device: Single-use, non-sterile medical mask designed to cover the nose and mouth of the healthcare professionals and/or the patient, including pediatric services, during surgical procedures, care or medical examinations, in order to limit the transmission of infectious agents and cross-contamination. Class I medical device.

Intended use :

Type II	Type IIR
Type II medical masks are intended to be worn by the healthcare professional during a patient examination, patient care or surgical procedure, in order to limit the transmission of infectious agents from medical staff to the patient. Medical masks can also be worn by patients or others to reduce the risk of spreading infection, particularly in an epidemic or pandemic situation.	Type IIR medical masks are intended to be worn by the healthcare professional during a patient examination, patient care or surgical procedure, in order to limit the transmission of infectious agents from medical staff to the patient and to protect the user of the mask against splashes of potentially contaminated biological fluids and viable particles present in the surgical environment. Medical masks can also be worn by patients or others to reduce the risk of spreading infection, particularly in an epidemic or pandemic situation.

Instructions for use : Before use, check the integrity of the mask, its links (earloops or tie-one) and any additional features such as visor, anti-fog, anti-glare or foam where applicable. Do not remove the mask from its packaging until ready for use. Wash your hands before use.

Fitting: Follow the instructions illustrated on the side of the box. The colour and/or marking of the mask is indicated on the outside of the device. The nose clip must be positioned upwards and shaped on the bridge of the nose to ensure a good seal. The mask must cover the nose, mouth and chin.

Warning: Do not use the product if its packaging is damaged. Discard and replace the mask:

- if it is damaged
- if it appears to interfere with breathing
- if it is contaminated by blood or any other infectious element or dirt.

If you experience any discomfort during use (irritation, warmth, nausea, dizziness, shortness of breath), remove the mask, immediately. If symptoms persist, contact your doctor. For use in a healthcare environment, dispose of the mask in the waste from healthcare activities involving infectious risks. Use on healthy skin only. Medical masks for children must be used under the supervision of an adult. Single-use device: do not reuse, do not wash the mask. Reusing a used device creates risks of contamination of the wearer or people in the vicinity due to the degradation of the mask's.

Limits of use: Medical masks are designed to be used for a maximum of 4 hours. Failure to comply with the instructions and warnings may reduce the effectiveness of the mask. This product does not eliminate the risk of contracting diseases or infections. This product is not personal protective equipment (PPE) meeting the requirements of Regulation (EU) 2016/425. This product contains a metal component (nose clip) and is not compatible for use in medical imaging.

Storage: Store the product in its original packaging, away from moisture, sunlight and light, in a cool environment or at room temperature.

Shelf life: The expiry date is indicated on the packaging.

Destruction: In accordance with current regulations.

Incident notification: Any serious incident occurring in connection with the use of the product must be notified to the manufacturer and to the competent authority of the Member State in which the user is established.

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