


1140
KOLMI COVERFEEL TOUCH - POWDER FREE NITRILE GLOVE

Class I Medical Device and Category III Personal Protective Equipment




FEATURES & BENEFITS

- Without accelerator to limit type IV allergies.
- Ultra comfortable glove and faster donning thanks to its internal polymer coating.
- Great elasticity and good dexterity.
- Non-slip texture for excellent grip and manipulation of objects.
- Waterproof and airtight glove.


INTENDED USE

- Glove intended for professional use.
- In a medical environment, this glove is used to provide a barrier against micro-organisms, between healthcare personnel, the patient, and their environment. Protection tested (permeation) against many chemotherapy agents (oncology).
- In an industrial environment, this glove is intended to protect the wearer from the tested chemicals (type B) and microorganisms (bacteria, molds and viruses).
- Powder-free, it reduces the risk of airborne contamination and limits the risk of irritation.
- Without natural latex, this glove is recommended for people allergic to it.
- This glove is suitable for food contact for any type of food.

TECHNICAL DESCRIPTION

Name of product:	Kolmi Coverfeel Touch
Type of product:	Single-use glove, no-sterile
Material:	Nitrile
Color:	Blue
Shape:	Ambidextrous
Edge :	Rolled
Texture :	Fingers textured
Powder :	No
Thickness of the palm:	0.07mm (±0,02mm)
Thickness of finger:	0.08mm (±0,02mm)
Length of glove:	≥ 240mm
	≥ 250 mm for XL & XXL
Unit weight (M):	3.5g
Country of manufacture:	FRANCE   
Conditioning :	Box from 100uc – Sizes XS, S, M & L
	Box from 90uc – Sizes XL & XXL

Article	Size	
1140-A	T6	XS
1140-B	T7	S
1140-C	T8	M
1140-D	T9	L
1140-E	T10	XL
1140-F	T11	XXL

CERTIFICATIONS & STANDARDS

- Designed and distributed under an ISO 13485 & ISO 9001 certified quality system.
- Complies with the requirements of Regulation (EU) 2017/745, relating to medical devices.
- Meets EN 455-1/2/3/4 standards for non-reusable medical gloves.
- Complies with the requirements of Regulation (EU) 2016/425, relating to personal protective equipment.
- Meets EN ISO 21420, EN ISO 374-1/2/4/5, EN 16523-1, ISO 16604 standards, for protection against tested chemicals (type B), microorganisms including viruses.
- In accordance with requirements of Reg.(EC) 1935/2004 and the French decree of 5th August 2020. This glove is suitable for food contact, with all types of foodstuff.
- Tested according to ASTM D6978-05 for the evaluation of the resistance of medical gloves to permeation by tested chemotherapy drugs.
- Electrostatic properties: tested to EN 1149-1, EN 1149-2 & EN 1149-3 (results on request).

Chemical permeation data (EN 16523-1:2015+A1:2018)		
Chemicals tested	N° CAS	Performance level/ Breakthrough time
J: n-Heptane	142-82-5	Level 3, < 120 min
K: Sodium hydroxide 40%	1310-73-2	Level 6, > 480 min
T: Formaldehyde 37%	50-00-0	Level 5, > 480 min

Resistance to cytotoxic drugs (ASTM D6978-05)	
Cytotoxic molecules and concentration	Minimum hole detection time
Carmustine (BCNU) - 3.3 mg/ml	21,7 min
Doxorubicin HCl - 2.0 mg/ml	> 240 min
Cyclophosphamide - 20 mg/ml	> 240 min
Paclitaxel - 6.0 mg/ml	> 240 min
Thiotepa - 10.0 mg/ml	54,3 min
Etoposide - 20.0 mg/ml	> 240 min
Fluorouracil - 50.0 mg/ml	> 240 min
Bulsulfan - 6mg/ml	> 240 min
Irinotecan HCl - 20mg/ml	> 240 min
Vidaza (Azacytidine) - 25 mg/ml	> 240 min
Fentanyl citrate injection - 100mcg/2ml	> 240 min
Methotrexate - 25mg/ml	> 240 min
Vincristine Sulfate - 1mg/ml	> 240 min
Ifosfamide - 50mg/ml	> 240 min
Cytarabine HCl - 100mg/ml	> 240 min



PRECAUTIONS

The device should be used on healthy skin only.

Single use, do not reuse.

After use, comply with the national regulations in force for the disposal of the device.

COMPOSITION OF THE DEVICE

Nitrile – without sulfide or accelerator

Latex: no

BPA: no

Phthalates CMR: no

Product of animal or plant origin: no

Presence of ethylene oxide: no

STORAGE CONDITIONS

Normal conservation and storage conditions: must not be exposed to humidity and sunlight.

Product lifespan: 3 years.

MARKING ON THE PACKAGING



LOGISTICS INFORMATION

CARTON SPECIFICATIONS				BOX SPECIFICATIONS			
Article	SIZE (mm)	Weight (kg)	Qty/pallet	BARCODE	SIZE (mm)	QTY / carton	BARCODE
1140-A	L344 x W256 x H270	3,2	60 cartons (6 layers of 10 cartons)	3 662 036 020 758	L250 x W121 x H64	10 dispenser boxes of 100 units	3 662 036 020 901
1140-B		3,7		3 662 036 020 765			3 662 036 020 918
1140-C		4,2		3 662 036 020 772			3 662 036 020 925
1140-D		4,7		3 662 036 020 789			3 662 036 020 932
1140-E		5		3 662 036 020 796		3 662 036 020 949	
1140-F		5,4		3 662 036 020 963		3 662 036 020 956	

