

1142

KOLMI COVERFEEL CARE - POWDER FREE NITRILE GLOVE

Class I Medical Device and Category III Personal Protective Equipment

FEATURES & BENEFITS

- Without accelerator to limit type IV allergies.
- Ultra comfortable glove and faster donning thanks to its internal polymer coating.
- Great elasticity and good dexterity.
- Robust for greater durability.
- Non-slip texture for excellent grip and manipulation of objects.
- Waterproof and airtight glove.



INTENDED USE

- Glove intended for professional use.
- In a medical environment, this glove is used to provide a barrier against micro-organisms, between healthcare personnel, the patient and their environment. Protection tested (permeation) against many chemotherapy agents (oncology).
- In an industrial environment, this glove is intended to protect the wearer from tested chemicals (type B), radioactive particles and microorganisms (bacteria, mold and viruses).
- Powder-free, it reduces the risk of airborne contamination and limits the risk of irritation.
- Without natural latex, this glove is recommended for people allergic to it.
- This glove is suitable for food contact for any type of food.

TECHNICAL DESCRIPTION

Name of product:	Kolmi Coverfeel Care
Type of product:	Single-use glove, no sterile
Material:	Nitrile
Color:	Blue
Shape:	Ambidextrous
Edge :	Rolled
Texture :	Fingers textured
Powder :	No
Thickness of the palm:	0.10mm (±0,02mm)
Thickness of finger:	0.11mm (±0,02mm)
Length of glove:	≥ 240mm ≥ 250 mm for XL & XXL
Unit weight (M):	5.2g
Country of manufacture:	FRANCE  
Conditioning :	Box from 100 uc –Sizes S, M & L Box from 90uc- Sizes XL & XXL

Article		Size	
1142-B	T7	S	
1142-C	T8	M	
1142-D	T9	L	
1142-E	T10	XL	
1142-F	T11	XXL	

CERTIFICATIONS & STANDARDS

- Designed and distributed under an ISO 13485 & ISO 9001 certified quality system.
- Complies with the requirements of Regulation (EU) 2017/745, relating to medical devices.
- Meets EN 455-1/2/3/4 standards for non-reusable medical gloves.
- Complies with the requirements of Regulation (EU) 2016/425, relating to personal protective equipment.
- Meets EN ISO 21420, EN ISO 374-1/2/4/5, EN 16523-1, ISO 16604 standards for chemicals tested (type B) and microorganisms including viruses. EN 421 for protection against radioactive particles.
- In accordance with requirements of Reg.(EC) 1935/2004 and the French decree of 5th August 2020. This glove is suitable for food contact, with all types of foodstuff.
- Tested according to ASTM D6978-05 for the evaluation of the resistance of medical gloves to permeation by the tested chemotherapy drugs.
- Electrostatic properties: tested to EN 1149-1, EN 1149-2 & EN 1149-3 (results on request).
- Compliant with ISO 18889:2019, relating to the protection of operators handling pesticides – Category G1. Test reports are available upon request.

Chemical permeation data (EN 16523-1:2015+A1:2018)

Chemicals tested	N° CAS	Performance level/ Breakthrough time
J: n-Heptane	142-82-5	Level 5, > 480 min
K: Sodium hydroxide 40%	1310-73-2	Level 5, > 480 min
T: Formaldehyde 37%	50-00-0	Level 5, > 480 min
Sulphuric acid 50%	7664-93-9	Level 5, > 480 min
Sodium hydroxide 50%	1310-73-2	Level 5, > 480 min
Peracetic acid 5%	79-21-0	Level 5, > 480 min
Cyclohexane 99%	100-82-7	Level 4, < 240 min
Isopropyl alcohol 70%	67-63-0	Level 3, < 120 min
P: Hydrogen peroxide 30%	7722-84-1	Level 2, < 60 min
Ethanol 70%	64-17-5	Level 2, < 60 min
Hydrochloric acid 37%	7647-01-0	Level 1, < 30 min
A: Methanol	67-56-1	Level 0, < 10 min
C: Acetonitrile	75-05-8	Level 0, < 10 min
L: Sulfuric Acid	7664-93-9	Level 0, < 10 min
M: Nitric Acid	7697-37-2	Level 0, < 10 min
N: Acetic Acid	64-19-7	Level 0, < 10 min
O: Ammonium hydroxide 25%	1336-21-6	Level 0, < 10 min
S: Hydrofluoric acid 40%	7664-39-3	Level 0, < 10 min
Methyl ethyl ketone (MEK) 100%	78-93-3	Level 0, < 10 min

Resistance to cytotoxic drugs (ASTM D6978-05)

Cytotoxic molecules and concentration	Minimum hole detection time
Cisplatin - 1mg/ml	> 240 min
Dacarbazine - 10mg/ml	186,6 min
Carmustine (BCNU) - 3.3 mg/ml	54,7 min
Doxorubicin HCl - 2.0 mg/ml	> 240 min
Cyclophosphamide - 20 mg/ml	> 240 min
Paclitaxel - 6.0 mg/ml	> 240 min
Thiotepa - 10.0 mg/ml	65,5 min
Etoposide - 20.0 mg/ml	> 240 min
Fluorouracil - 50.0 mg/ml	> 240 min
Bulsulfan - 6mg/ml	> 240 min
Irinotecan HCl - 20mg/ml	> 240 min
Vidaza (Azacytidine) - 25 mg/ml	> 240 min
Fentanyl citrate injection - 100mcg/2ml	> 240 min
Methotrexate - 25mg/ml	> 240 min
Vincristine Sulfate - 1mg/ml	> 240 min
Ifosfamide - 50mg/ml	> 240 min
Cytarabine HCl - 100mg/ml	> 240 min

PRECAUTIONS

The device should be used on healthy skin only.
Single use, do not reuse.
After use, comply with the national regulations in force for the disposal of the device.


COMPOSITION OF THE DEVICE

Nitrile – without accelerator
Latex: no
BPA: no
Phthalates CMR: no
Product of animal or plant origin: no
Presence of ethylene oxide: no

STORAGE CONDITIONS

Normal conservation and storage conditions: must not be exposed to humidity and direct sunlight.
Product lifespan: 3 years.

MARKING ON THE PACKAGING



LOGISTICS INFORMATION

CARTON SPECIFICATIONS					BOX SPECIFICATIONS		
Article	SIZE (mm)	Weight (kg)	Qty/pallet	BARCODE	SIZE (mm)	QTY / carton	BARCODE
1142-B	L390 x l290 x H275	5,7	48 boxes	3 662 036 021 045	L255 x W135 x H73	10 dispenser boxes of 100 units	3 662 036 020 987
1142-C		6		3 662 036 021 052			3 662 036 020 994
1142-D		6,5		3 662 036 021 069			3 662 036 021 007
1142-E		6,7		3 662 036 021 076		10 dispenser boxes of 90 units	3 662 036 021 014
1142-F		7,1		3 662 036 021 083			3 662 036 021 021