

10.007/XXL

**WHITE PP SHOECOVERS
WITH WHITE PLASTIC COATED SOLE – XXL SIZE**

Medical device Class I

FEATURES & BENEFITS

- PP shoecovers with 60µm PE plastic-coated sole consist of a PE non-slip preventing users from slipping and ensures a long product life longer face the risk of abrasion on rough floors
- Available for color code tracking
- Adaptable to different morphologies of all types of shoes
- Due to its large size and height, this model fits perfectly to all types of safety shoes



PRODUCT DESCRIPTION

Product name : White PP shoecovers SafeFeet® SkidGuard™ with white plastic-coated sole

Product type : Single-use, non sterile

Material : PP spunlaid nonwoven + PE

Confection : Gluing

Grammage : 35 g/m²

Sole thickness : 60µm

Size : 230 x 420 mm – XXL Size & very tall height

Unit weight : 14,6 g ■ ■

Country of origin : France

INDUSTRIES

Industry / Chemical / / Food / Away from home food / Non residential catering / Community / Administration / Petrochemical Health / Medical / Surgical / Sanitary / Maintenance / Cleaning / Clean room / Lab / Pharma / Research

MECHANICAL INFORMATION					
Measurements done by independant laboratories (measures are available upon request)					
TEST	STANDARD		AVERAGE	UNIT	TOLERANCE
Grammage	ISO 9073-1		35	g/m ²	±10%
Tensilestrenght	ISO 9073-3	Tensilestren. MD:	64,90	N	±10%
		Tensilestren. CD:	80,20	N	±10%
Elongation		Elongation MD:	50,74	%	±10%
		Elongation CD:	63,24	%	±10%

10.007/XXL

WHITE PP SHOECOVERS WITH WHITE PLASTIC COATED SOLE – XXL SIZE

Medical device Class I

MICROBIOLOGICAL INFORMATION

Assessment of initial microbial contamination according to standards NF EN ISO 11737:2018 + A1: 2021.
Additional microbiological controls: ASR, E.coli, Staphylococci.

CERTIFICATIONS & STANDARDS

Meet the requirements of Regulation (EU) 2017/745 on medical devices.
Complies with applicable standards EN ISO 14971.
Manufacturing site : ISO 9001 – ISO 13485.

PRECAUTIONS

The product must be used on intact skin only.
The reuse or extended use may lead to infection or cross-contamination.
Please refer to your national legislation on the subject of disposal of this product once it has been used.

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.
Shelf life of the product: 5 years.

MARKING ON PACKAGING



* Do not use if package is damaged.

INFORMATIONS LOGISTIQUES

OUTER CASE SPECIFICATIONS					BAG SPECIFICATIONS
Article	SIZE (mm)	WEIGHT (kg)	Qty / pallet	BARCODE	Packaging
10.007/XXL	400 x 400 x 375	5,0	36 cartons	3 700 381 503 264	1 bag of 300 pcs