

1125 LATEX GLOVE SAFETOUCH CONNECT POWDERED

Personal Protective Equipment, Category II

FEATURES & BENEFITS

- Great elasticity
- Agreeable dexterity
- Powdered, its donning is easier
- Watertight and airtight
- Contraindication: Latex is a natural material that contains protein which might cause allergies



Ideal for:

- Glove for professional use
- In the healthcare and industrial environment, this glove is used to make barrier against microorganisms between the wearer of the glove and the examined people or objects / materials handled.
- For food handling / preparations (HACCP)

PRODUCT DESCRIPTION

Product name :	SafeTouch Connect
Product type :	Single use, non-sterile
Material :	Latex
Colour :	Natural
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Smooth
Powdered :	Yes
Length :	240 mm
Unit Weight (M):	6±0,2 g
Palm thickness :	0,10±0,02mm
Finger thickness :	0,12±0,02mm
Country of origin :	Malaysia

Item	Size
1125-A	T 6 XS
1125-B	T 7 S
1125-C	T 8 M
1125-D	T 9 L
1125-E	T 10 XL

Marking on packaging:



* Do not use if package is damaged

INDUSTRIES

Laboratory / Community / Administration / Health / Medical / Sanitary / Cleaning / Pharmaceutical / Research.

CERTIFICATION

Manufactured under ISO 13485 & ISO 9001 certified quality system.

In accordance with requirements of European Reg. 2016/425 (Personal protective Equipments).

Comply with the EN 420 & EN374-2/5(*) for protection against micro-organisms.

In accordance with requirements of Reg.(EC) 1935/2004 and French decree of 9th November, 1994. This glove is suitable for food contact**.

* Uncontrolled against viruses.

**Restriction : not put in contact with acidic foods (pH<4,5).

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be exposed to moisture and sunlight, must be stored at a temperature between 5°C and 40°C.

Shelf life of the product: 5 years

LOGISTIC INFORMATION








Specifications outer case



Specifications inner case

Colour code dispenser and carton

Item	Size mm	Gross weight kg	QTY/pallet	Size mm	QTY	Colour
1125-A	360 x 250 x 230	5,86	63	215 x 120 x 70	10 x 100	
1125-B	360 x 250 x 230	6,51	63	215 x 120 x 70	10 x 100	
1125-C	360 x 250 x 230	6,92	63	215 x 120 x 70	10 x 100	
1125-D	360 x 250 x 230	7,23	63	215 x 120 x 70	10 x 100	
1125-E	360 x 250 x 230	7,23	63	215 x 120 x 70	10 x 100	

EAN Codes

Item	Inner case	Outer case
1125-A	3 700 381 514 314	3 700 381 514 260
1125-B	3 700 381 514 321	3 700 381 514 277
1125-C	3 700 381 514 338	3 700 381 514 284
1125-D	3 700 381 514 345	3 700 381 514 291
1125-E	3 700 381 514 352	3 700 381 514 307