

# 1163

**LATEX GLOVE SAFETOUCH CONNECT REJUVENATE POWDER-FREE**  
**Enriched in Lanolin and Vitamin E**  
**Class I, Medical Device**

## FEATURES & BENEFITS

- Great elasticity
- Good dexterity
- Latex is the material which allows the best contact sensibility
- Its micro-textured appearance gives it correct touch and grip of objects
- Enriched in Lanolin and Vitamin E, this glove gives the feeling of a better hydrated skin, softer and smoother.
- Contraindication: Latex is a natural material that contains protein which might cause allergies



## Recommendation:

- In the healthcare environment, this glove is used to make barrier against microorganisms between the wearer of the glove and the examined people or objects / materials handled.
- Powder-free, it reduces cross contamination and limits the risk of irritation.

## PRODUCT DESCRIPTION

Product name :	<b>SafeTouch ADVANCED Rejuvenate</b>
Product type :	Single use, non-sterile
Material :	Latex
Colour :	Natural
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Micro-textured
Powdered :	No
Length :	240 mm
Palm thickness :	0,11±0,03 mm
Finger thickness :	0,13±0,03 mm
Integrated lotion:	Lanolin and Vitamin E
Unit weight (M):	5,2g (±0,2g)
Country of origin :	Malaysia

Item	Size
1163-A	T 6 XS
1163-B	T 7 S
1163-C	T 8 M
1163-D	T 9 L
1163-E	T10 XL

Marking on packaging:



\* Do not use if package is damaged

## INDUSTRIES

Health / Medical

## CERTIFICATION

Manufactured under an ISO 13485 and ISO 9001 certified quality system.  
Meets the requirements of European Regulation 2017/745 (Medical Devices).  
Comply with the standards 455-1/2/3/4\_Medical gloves for single use.

## PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.  
The reuse or extended use may lead to infection or cross-contamination.  
Please refer to your national legislation on the subject of disposal of this product once it has been used.

## STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C.  
Shelf life of the product: 5 years.

## LOGISTIC INFORMATION



### Specifications outer case

Item	Size mm	Gross weight kg	QTY/ Pallet
1163-A	370 x 255 x 250	6,3	54
1163-B	370 x 255 x 250	6,8	54
1163-C	370 x 255 x 250	7,1	54
1163-D	370 x 255 x 250	8,3	54
1163-E	370 x 255 x 250	8,8	54



### Specifications inner case

Size mm	QTY
240 x 120 x 70	10 x 100
240 x 120 x 70	10 x 100
240 x 120 x 70	10 x 100
240 x 120 x 70	10 x 100
240 x 120 x 70	10 x 100

### Colour code inner and outer

Colour






### EAN Codes

Item	Inner case	Outer case
1163-A	3 700 381 514 970	3 700 381 514 581
1163-B	3 700 381 514 987	3 700 381 514 598
1163-C	3 700 381 514 994	3 700 381 514 604
1163-D	3 700 381 515 007	3 700 381 514 611
1163-E	3 700 381 520 285	3 700 381 520 292