

1174N SAFETOUCH ADVANCED PLATINUM POWDER-FREE NITRILE GLOVE

Class I Medical Device and Category III Personal Protective Equipment

FEATURES & BENEFITS

- Great elasticity.
- Good dexterity.
- Its smooth appearance gives it correct touch and grip of objects.
- Ultra comfortable, this generation allows an excellent stretchability and an optimal comfort.
- Watertight and airtight.



RECOMMENDATION:

- Glove for professional use
- In the healthcare environment, this glove is used to make barrier against microorganisms between the professional and the patient and their environment.
- In an industrial environment, this glove is intended to protect the wearer from tested chemicals and microorganisms (bacteria, moulds and viruses).
- The powder-free glove reduces the risk of airborne contamination and limits the risk of irritation.

PRODUCT DESCRIPTION

Product name :	SafeTouch Advanced Platinum
Product type :	Single use, non-sterile
Material :	Nitrile
Colour :	White
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Textured fingertips
Powdered :	No
Palm thickness :	0,07 mm (± 0.02 mm)
Finger thickness :	0,10 mm (± 0.02 mm)
Unit weight (M):	3,5 g (± 0,2g)
Country of origin :	Malaysia

Article	Size	
1174N-A	T 6	XS
1174N-B	T 7	S
1174N-C	T 8	M
1174N-D	T 9	L
1174N-E	T 10	XL

Marking on packaging:



INDUSTRIES

Health / Medical / Dental / Industry/ Food

1174N SAFETOUCH ADVANCED PLATINIUM POWDER-FREE NITRILE GLOVE

Class I Medical Device and Category III Personal Protective Equipment

CERTIFICATION

Manufactured under an ISO 13485 and ISO 9001 certified quality system.

In accordance with requirements of Reg (EU) 2017/745, UK Medical Device Regulation 2002 and the Swiss Ordonnance (ODim) from July 1st 2020, on medical devices.

Comply with the standards EN 455-1/2/3/4 for single use medical gloves for single use.

In accordance with requirements of Reg (EU) 2016/425 and Regulation 2016/425 as amended to apply in GB, on personal protective equipments.

Comply with the EN 420, EN ISO 374-1/2/4/5, for protective gloves against tested chemicals and microorganisms including viruses.

In accordance with requirements of Reg.(EC) 1935/2004 and the French decree of 5th August, 2020. This glove is suitable for food contact, with all types of foodstuff.

PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.

The reuse or extended use may lead to infection or cross-contamination.

Please refer to your national legislation on the subject of disposal of this product once it has been used.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, store in a cool environment or at room temperature.

Shelf life of the product: 3 years.

LOGISTIC INFORMATION








Specifications outer case



Specifications inner case

Colour code inner and outer

Item	Size mm	Gross weight kg	QTY/ pallet	Size mm	QTY	Colour
1174N-A		4,1				
1174N-B		4,7				
1174N-C	262 x 252 x 270	5,1	60 cartons (6 layers of 10 cartons)	245 x 125 x 48	10 boxes of 100 pieces	
1174N-D		5,3				
1174N-E		6,2				



EAN Codes

Item	Outer case	Inner case
1174N-A	3 700 381 519 661	3 700 381 519 678
1174N-B	3 700 381 519 685	3 700 381 519 692
1174N-C	3 700 381 519 708	3 700 381 519 715
1174N-D	3 700 381 519 722	3 700 381 519 739
1174N-E	3 700 381 519 746	3 700 381 519 753