

1176 NITRILE GLOVE SAFETOUCH ADVANCED REJUVENATE POWDER-FREE

Enriched in Lanolin and Vitamin E

Class I, Medical Device

FEATURES & BENEFITS

- Great elasticity
- Good dexterity
- Its micro-textured appearance gives it correct touch and grip of objects
- Enriched in Lanolin and Vitamin E, this glove gives the feeling of a better hydrated skin, softer and smoother.



Recommendation:

- In the healthcare environment, this glove is used to make barrier against microorganisms between the wearer of the glove and the examined people or objects / materials handled.
- Powder-free, it reduces cross contamination and limits the risk of irritation.

PRODUCT DESCRIPTION

Product name :	SafeTouch ADVANCED Rejuvenate - Glove
Product type :	Single use, non-sterile
Material :	Nitrile
Colour :	White
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Micro-textured
Powdered :	No
Length :	240 mm
Palm thickness :	0,09±0,03 mm
Finger thickness :	0,11±0,03 mm
Integrated lotion:	Lanolin and Vitamin E
Unit weight (M):	4g (±0,2g)
0,2g)Country of origin :	Malaysia

Item	Size
1176-A	T 6 XS
1176-B	T 7 S
1176-C	T 8 M
1176-D	T 9 L
1176-E	T10 XL

Marking on packaging:        *

*Do not use if package is damaged

INDUSTRIES

Health / Medical

CERTIFICATION

Manufactured under an ISO 13485 and ISO 9001 certified quality system.
Meets the requirements of European Regulation 2017/745 (Medical Devices).
Comply with the standards 455-1/2/3/4_Medical gloves for single use.








PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.
The reuse or extended use may lead to infection or cross-contamination.
Please refer to your national legislation on the subject of disposal of this product once it has been used.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C.
Shelf life of the product: 5 years.

LOGISTIC INFORMATION

 Specifications outer case				 Specifications inner case		Colour code inner and outer
Item	Size mm	Gross weight kg	QTY/ Pallet	Size mm	QTY	Colour
1176-A	320 x 265 x 250	4,4	63	235 x 130 x 60	10 x 100	
1176-B	320 x 265 x 250	4,7	63	235 x 130 x 60	10 x 100	
1176-C	320 x 265 x 250	5,1	63	235 x 130 x 60	10 x 100	
1176-D	320 x 265 x 250	5,3	63	235 x 130 x 60	10 x 100	
1176-E	320 x 265 x 250	5,5	63	235 x 130 x 60	10 x 100	

EAN Codes

Item	Inner case	Outer case
1176-A	3 700 381 514 802	3 700 381 514 413
1176-B	3 700 381 514 819	3 700 381 514 420
1176-C	3 700 381 514 826	3 700 381 514 437
1176-D	3 700 381 514 833	3 700 381 514 444
1176-E	3 700 381 521 565	3 700 381 521 558