

## 1193N NITRILE GLOVE SAFETOUCH ADVANCED VITALS POWDER-FREE Class I Medical Device and Category III Personal Protective Equipment

### FEATURES & BENEFITS

- Latex-like feel.
- Agreeable dexterity.
- Excellent grip in both wet and dry conditions.
- Watertight and airtight.



### RECOMMENDATION:


- In the healthcare and industrial environment, this glove is used to make barrier against microorganisms between the wearer of the glove and the examined people or objects/ materials handled.
- Powder-free, it reduces cross contamination and limits the risk of irritation.
- In an industrial environment, this glove is used to protect the wearer from microorganisms (bacteria, molds and viruses) and tested chemicals.
- For food handling / preparations (HACCP).

### PRODUCT DESCRIPTION

Product name:	<b>SafeTouch® Advanced</b>
Product type :	Single use, non-sterile
Material :	Nitrile
Colour :	Blue
Shape :	Ambidextrous
Cuff:	Rolled
Surface :	Textured fingertips
Powdered :	No
Length :	240 mm
Palm thickness :	0,07 mm (±0,02mm)
Finger thickness :	0,09 mm (±0,02mm)
Unit weight (M):	3,2 g
Country of origin :	Malaysia

Item	Size
1193N-A	T 6 XS
1193N-B	T 7 S
1193N-C	T 8 M
1193N-D	T 9 L
1193N-E	T 10 XL

Marking on packaging:



### INDUSTRIES

Food / Health / Medical / Dental/ Administration / Hygiene/ Industry

## CERTIFICATION

Manufactured under an ISO 13485 and ISO 9001 certified quality system.

In accordance with requirements of Reg (EU) 2017/745 and the Swiss Ordonnance (ODim) from July 1st 2020, on medical devices.

Comply with the standards EN 455-1/2/3/4 for single use medical gloves.

In accordance with requirements of Reg (EU) 2016/425, on personal protective equipments.

Comply with the EN ISO 21420, EN ISO 374-1/2/4/5, for protective gloves against tested chemicals and microorganisms including viruses.

In accordance with requirements of Reg.(EC) 1935/2004 and the French decree of 5th August, 2020. This glove is suitable for food contact, with all types of foodstuff (see annex 1).

## PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.

The reuse or extended use may lead to infection or cross-contamination.

Please refer to your national legislation on the subject of disposal of this product once it has been used.

## STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored in a cool environment or at room temperature.

Shelf life of the product: 3 years.

## LOGISTIC INFORMATION








### Specifications outer case



### Specifications inner case

### Colour code inner and outer

Item	Size mm	Gross weight kg	QTY/pallet	Size mm	QTY	Colour
1193N-A	300 x 260 x 235	3,2	60	245 x 125 x 48	10 x 100	
1193N-B	300 x 260 x 235	3,4	60	245 x 125 x 48	10 x 100	
1193N-C	300 x 260 x 235	3,7	60	245 x 125 x 48	10 x 100	
1193N-D	300 x 260 x 235	4,5	60	245 x 125 x 48	10 x 100	
1193N-E	300 x 260 x 235	5,1	60	245 x 125 x 48	10 x 100	



### EAN Codes

Item	Inner case	Outer case
1193N-A	3 700 381 519 999	3 700 381 519 982
1193N-B	3 700 381 520 018	3 700 381 520 001
1193N-C	3 700 381 520 032	3 700 381 520 025
1193N-D	3 700 381 520 056	3 700 381 520 049
1193N-E	3 700 381 520 070	3 700 381 520 063