

30.350N NITRILE SAFETOUCH ADVANCED X' TRA BLUE GLOVE – POWDER FREE


Personal Protective Equipment, Category III

FEATURES & BENEFITS

- Feeling like latex
- Good dexterity
- Long oversleeve 355 mm
- Textured fingertips to assure a better grip in wet and dry conditions
- Provide protection against micro-organisms, radioactive contamination and protects the wearer from tested chemicals
- Internal /External chlorination to assist threading

Recommendation:

- Users allergic to latex
- Contact with chemical products (cf instruction for use)
- Nuclear laboratories: pre-glove application for glove box



PRODUCT DESCRIPTION

Product name :	SafeTouch Advanced X'tra
Product type :	Single Use non-sterile
Material :	Nitrile
Colour :	Blue
Shape :	Ambidextrous
Edges :	Rolled
Surface :	Finger textured
Powder:	Powder free
Length (mm) :	355 mm
Palm thickness:	0,16 mm (+/- 0,02)
Finger thickness :	0,21 mm (+/- 0,02)
Weight (M) :	13,1 g
Country of origin:	Malaysia

Item	Size	
30.350N/S	T 7	S
30.350N/M	T 8	M
30.350N/L	T 9	L
30.350N/XL	T 10	XL

Marking on packaging:



INDUSTRIES:

Industry/ Chemical / Petrochemical / Health / Medical / Maintenance / Cleaning / Clean room / Laboratory / Pharma / Research

30.350N

NITRILE SAFETOUCH ADVANCED X'TRA BLUE GLOVE – POWDER FREE

Personal Protective Equipment, Category III

STANDARDS AND CERTIFICATION

Manufactured under an ISO 13485 and ISO 9001 certified quality system.
 Meets the requirements of European Regulation 2016/425 (Personal Protective Equipments).
 Meets EN ISO 21420 : 2020 - General requirements for protective gloves.
 Meets EN 421: 2010 - Requirements for protective gloves against radioactive contamination.
 Doesn't protect against ionizing radiation and is used outside containment.
 Meets EN 374-1 / 2/4/5*, EN 16523-1 standards for protective gloves against chemicals tested and micro-organisms.
 *Non tested against viruses.

PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.
 The reuse or extended use may lead to infection or cross-contamination.
 Please refer to your national legislation on the subject of disposal of this product once it has been used.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C.
 Shelf life of the product: 5 years.

LOGISTIC INFORMATIONS



Specification carton (double groove)

Specification polybag

Color code polybag & carton

Item	Size mm	Gross weight	Carton / pallet	Size mm	Packaging	Color
30.350/S	140 x 430 x 260	6	56	430 x 310 x 20	5 x 100	
30.350/M	140 x 430 x 260	6,9	56	430 x 310 x 20	5 x 100	
30.350/L	140 x 430 x 260	7	56	430 x 310 x 20	5 x 100	
30.350/XL	140 x 430 x 260	7,4	56	430 x 310 x 20	5 x 100	



EAN code

Item	Carton	Polybag	Item	Carton	Polybag
30.350/S	3 700 381 520 322	3 700 381 520 339	30.350/L	3 700 381 520 360	3 700 381 520 377
30.350/M	3 700 381 520 346	3 700 381 520 353	30.350/XL	3 700 381 520 384	3 700 381 520 391