

ORANGE PE GLOVE, TRANSLUCENT 25µm, 920 mm

PRODUCT DETAILS :



Product name : Orange PE Glove, Translucent 25µm, 920 mm
Material : polyethylene
AQL : 4
Colour : orange
Sterilization : no
Thickness: 25µm (fingers and sleeves)
Shape : Ambidextrous
Powdered : No
Lenght : 920 mm
Unit weight : 9 g
Origin : Asia



FEATURES & BENEFITS:

Design :

- Polyethylene offers a great dexterity
- It can be used wboth by the left hand and the right hand

LOGISTICS INFORMATION :



Packaging : Quantity : 100 pcs / polybag
Type of packaging : carton
Number of polybag / carton : 20 polybags
Quantity per carton : 2000 pcs

Carton :

Carton size (mm) : 350x250x360
Gross weight : 16,860 kg
Net weight : 16,560 kg

Polybag:

Polybag size (mm) : 210x190
Gross weight : 828 gr
Net weight : 725 gr

Palettization : 20 cartons of 20 polybags : 40 000 pcs
Pallet weight : 344 kg

CONDITIONS AND STORAGE CONDITIONS:

Normal conditions of conservation & storage: not to be exposed to moisture and sun

Special precautions: single use

Presence of indicators of temperature if it is necessary