

**50.401**  
**50.401BL**

**APRON PE 30 µm**  
WHITE or BLUE – 120 x 70 m – Dispenser bag wall

## FEATURES & BENEFITS

- Resistant PE tie bonds.
- Distribution system « 1 by 1 » to avoid the cross contamination of the entire batch when leaving aprons of their dispenser bag.
- Possibility of identification by color codes.
- Thanks to its large size, the apron offers better protection than a standard one because it provides better cover.
- Food contact.

## TECHNICAL DESCRIPTION

Product name:	Apron PE 30 µm
Product type :	Single use, non sterile
Material :	Polypropylene
Confection :	Cut / Folding
Color :	White or blue
Grammage/Thickness :	30 µm
Dimensions :	120 x 70 cm
Weight :	15 g
Packaging :	Dispenser bag wall
Country of origin :	China



**50.401**

**50.401BL**

## INDUSTRIES

Industrial / Restaurant & catering / Community services / Hygiene

## STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

## CERTIFICATIONS & STANDARDS

Manufacturing site : ISO 9001.

Meets the requirements of Reg.(EC)1935/2004 and Reg. (EU)10/2011. This product is suitable for food contact with all food types.

## INFORMATIONS LOGISTIQUES

OUTERCASE SPECIFICATIONS					SPECIFICATIONS CDT	
Article	SIZE (mm)	WEIGHT (kg)	Qty / pallet	EAN CODE	CDT	CODE EAN
50.401	475 x 375 x 230	13	20 outercases	3 700 381 500 973	10 polybags of 100 units	3 700 381 502 427
50.401BL				3 700 381 520 582		3 700 381 520 599