

813600

Surgical Caps

813700

Medical Device Class I

FEATURES & BENEFITS

- Airy and extensible non-woven
- Very aesthetic and very covering cap
- Instructions for use:
 1. To be more aesthetic, turn the seams
 2. Adjust the height by turning the edge

Recommendation :



- Can be used by nurses, surgeons or patients in the operating room



813700

TECHNICAL DESCRIPTION

Product name: MEDICOM Surgical caps – Surgical caps
 Product type: Single use, non sterile
 Material : 84% Viscose + 16% Ethylene copolymer– vinyl acetate
 Confection : Sewn
 Color : Blue / Green
 Dimensions : Circumference : 56 cm (± 2 cm)
 Height : 32 cm (± 2 cm)
 Unit weight : 7,2 g (± 10%)
 Packaging : Dispenser box of 100 units
 Origin : France

Item	Color
813600	Blue 
813700	Green 

Marking on packaging:



Mechanical Information (Measurements done by Cetelor Laboratory)

TEST	Standard		Average	Unit	Tolerance
Grammage	ISO 9073-1		40	g/m ²	± 10%
Tensile strength	ISO 9073-3	Dry MD:	55,4	N	± 10%
		Dry CD:	8,6	N	± 10%
Elongation	ISO 9073-3	Dry MD:	10,5	%	± 10%
		Dry CD:	71,1	%	± 10%

Manufacturing Microbiological Information (by ICARE Laboratory)

Bioburden evaluation meets the requirements of EN ISO 11737: 2018+A1 2021
 Additional microbiological controls: ASR, E.coli, Staphylococcus.

813600
813700

Surgical Caps
Medical Device Class I

STANDARDS & CERTIFICATION

Manufactured under ISO 13485 certified quality system.
In conformity with the relevant harmonized standard EN ISO 14971.
Complies with the requirements of Regulation (EU) 2017/745 relating to Medical Devices.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, store in a cool environment or at room temperature.
Product life: 5 years.

PRECAUTIONS

The device should only be used on clean, intact skin.
Their re-use or prolonged use can produce infection or cross-contamination.
After use, comply with the national regulations in force for the disposal of the device.

LOGISTICAL INFORMATION



Outer case specifications

Dispenser box specifications

Item	Size mm	Gross weight kg	QTY/pallet	Size mm	QTY
813600 813700	575 x 225 x 370	4,6	35 cases (5 layers of 7 cases)	200 x 108 x 350	5 dispenser boxes of 100 pieces



EAN Code

Item	Outer case	Dispenser box
813600	3 662 036 015 198	3 662 036 015 204
813700	3 662 036 015 211	3 662 036 015 228

