

## A69955

### WATERPROOF PP ISOLATION GOWN

With PE coating

#### FEATURES & BENEFITS

- Premium protection to meet all your infection control needs
- Comfort and resistant
- Large opening for the head
- Patented « Flexneck » technology with a flexible elastic design
- Long sleeves
- Thumbloop for extra protection

#### Recommended use:

- Disposable Impervious Isolation Gown designed to protect health professionals and patients from cross infection



#### PRODUCT DETAILS

|               |   |
|---------------|---|
| Product name: | Impervious coated PP Isolation Gown   |
| Product type: | Single use, non-sterile   |
| Material:     | Polypropylene + Polyethylene  |
| Confection:   | Ultrasonically sealed, without collar, without pocket, Flexneck neck closure, long sleeves + thumb loops, elasticated wrists, Back closure with laces |
| Colour:       | Blue  |
| Grammage:     | 40 g/m <sup>2</sup> (± 10%)   |
| Unit weight:  | 118 g (± 10%)   |
| Packaging:    | Polybag of 10 units   |
| Origin:       | China   |

#### INDUSTRIES

Health / Medical

| Item   | Size | Length<br>(± 1 cm) | Underarm width<br>(± 1 cm) | Gown width with sleeves<br>(includes thumb loop)<br>(± 1 cm) | Open gown width<br>(± 1 cm) |
|--------|------|--------------------|----------------------------|--|-----------------------------|
| A69955 | XL   | 129 cm             | 79 cm                      | 235 cm   | 158 cm                      |

#### CERTIFICATION & STANDARDS

In accordance with ANSI/AAMI PB 70 level II barrier performance for surgical gown and tested for chemotherapy drug permeation as per ASTM D6978.

## STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

## LOGISTICAL INFORMATION

| OUTERCASE SPECIFICATION |                 |             |  |                    | CDT SPECIFICATIONS      |                    |
|-------------------------|-----------------|-------------|--|--------------------|-------------------------|--------------------|
| Article                 | SIZE (mm)       | WEIGHT (kg) | Qty / pallet                             | EAN CODE           | CDT                     | CODE EAN           |
| A69955                  | 540 x 325 x 400 | 13          | 20 outercases (5 layers of 4 outercases) | 206 868 640 518 97 | 10 polybags of 10 units | 106 868 640 518 90 |