

# A69962 SMS ISOLATION GOWN

## FEATURES & BENEFITS

- Premium protection to meet all your infection control needs
- Comprised of one layer of meltblown polypropylene bonded between two layers of spunbond polypropylene
- 3-ply protective fabric that combines strength, fluid resistance with soft, breathable comfort
- Large opening for the head
- Patented « Flexneck » technology with a flexible elastic design
- Long sleeves
- Thumbloop for extra protection

### Recommended use:

- Disposable SMS Polypropylene Isolation Gown designed to protect health professionals and patients from cross infection
- To wear over clothing



## PRODUCT DETAILS

Product name:	SMS Isolation Gown
Product type:	Single use, non-sterile
Material:	SMS Polypropylene
Confection:	Ultrasonically sealed, without collar, without pocket, Flexneck neck closure, long sleeves + thumb loops, elasticated wrists, Back closure with laces
Colour:	Yellow
Grammage:	> 23,5 g/m <sup>2</sup>
Unit weight:	70 g (± 10%)
Packaging:	Polybag of 10 units
Origin:	China

## INDUSTRIES

Health / Medical

Item	Size	Length (± 1 cm)	Underarm width (± 1 cm)	Gown width with sleeves (includes thumb loop) (± 1 cm)	Open gown width (± 1 cm)
A69962	XL	129 cm	79 cm	209 cm	158 cm

## CERTIFICATION & STANDARDS

In accordance with ANSI/AAMI PB 70 level II barrier performance for surgical gown.

## STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

## LOGISTICAL INFORMATION

OUTERCASE SPECIFICATION					CDT SPECIFICATIONS	
Article	SIZE (mm)	WEIGHT (kg)	Qty / pallet	EAN CODE	CDT	CODE EAN
A69962	630 x 290 x 335	8	24 outercases (6 layers of 4 outercases)	206 868 640 519 58	10 polybags of 10 units	106 868 640 519 51