

FICHE TECHNIQUE

WIPES

Embossed Airlaid 50 g/m²

28 x 38 cm

Kolmi®



Item	Colour
E01303	<input type="checkbox"/> White


FEATURES & BENEFITS

- Very absorbent and resistant nonwoven sheet

PRECONISATION :

- Used for wiping body
- Used for surface cleaning

TECHNICAL DESCRIPTION

Product name: Wipes – Embossed Airlaid 50 g/m²
Product type: Single use, non sterile
Material : Embossed (2 sides) Airlaid nonwoven - cellulosic fibers
Grammage : 50 g/m²
Unit Weight: 5,8 g (±10%)
Size: 280 x 380 mm (±5mm)
Packaging: Bag of 50 units
Colour : White
Origin : France 

MARKING ON PACKAGING



* Do not use if package is damaged

MECHANICAL INFORMATION

TEST	TEST METHOD	AVERAGE	UNIT	TOLERANCE	
Grammage	EDANA 40.3-90	50	g/m ²	± 2g/m ²	
Thickness	EDANA 30.5-99	0,76	mm	± 0,07mm	
Tensile strength	EDANA 20.2.89	Dry MD :	13	N/5cm	± 1,5
		Dry CD :	11	N/5cm	± 1,1
		Wet MD :	8	N/5cm	± 0,8
		Wet CD :	5,5	N/5cm	± 0,5
Absorption	EDANA 10.4-02	12	g/g	± 10%	

MICROBIOLOGICAL INFORMATION

Assessment of initial microbial contamination in accordance with NF EN ISO 11737 : 2018 + A1 :2021.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTICAL INFORMATION

OUTER CASE SPECIFICATIONS					
Article	Size (mm)	Weight (kg)	Qty/ palet	Qty/ box	EAN code bag
E01303	600 x 400 x 219	4,7	40 boxes (4 boxes x 10 layers)	12 bags of 60 units	3 662 036 001 122