

FICHE TECHNIQUE

WIPES

Plain Airlaid 60 g/m²

32 x 33 cm

Kolmi®



Item	Colour
E02307	<input type="checkbox"/> White


FEATURES & BENEFITS

- Very absorbent and resistant nonwoven sheet

PRECONISATION :

- Used for wiping body
- Used for surface cleaning

TECHNICAL DESCRIPTION

Product name: Wipes – Plain Airlaid 60 g/m²
Product type: Single use, non sterile
Material : Plain Airlaid nonwoven - cellulosic fibers
Grammage : 60 g/m²
Unit Weight: 5,4 g (±10%)
Size: 320 x 330 mm (±5mm)
Packaging: Bag of 50 units
Colour : White
Origin : France 

MARKING ON PACKAGING



* Do not use if package is damaged

MECHANICAL INFORMATION

TEST	STANDARD		AVERAGE	UNIT	TOLERANCE
Grammage	ISO 9073-1		60	g/m ²	± 5%
Thickness	ISO 9073-2		0,76	mm	± 10%
Tensile strength	ISO 9073-3	Dry MD :	16	N/5cm	± 10%
		Dry CD :	11	N/5cm	± 10%
		Wet MD :	8	N/5cm	± 10%
		Wet CD :	5,5	N/5cm	± 10%
Absorption	ISO 9073-6		1200	%	± 10%

MICROBIOLOGICAL INFORMATION

Assessment of initial microbial contamination in accordance with NF EN ISO 11737 : 2018 + A1 :2021.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTICAL INFORMATION

OUTER CASE SPECIFICATIONS					
Article	Size (mm)	Weight (kg)	Qty/ palet	Qty/ box	EAN code bag
E02307	600 x 400 x 250	4	36 boxes (4 boxes x 9 layers)	12 bags of 50 units	3 662 036 001 221