

# FICHE TECHNIQUE

## WIPES

Embossed Airlaid 60 g/m<sup>2</sup>

28 x 38 cm / 28 x 30 cm / 49 x 45 cm

# Kolmi®




### FEATURES & BENEFITS

- Very absorbent and resistant nonwoven sheet

### PRECONISATION:

- Used for wiping body
- Used for surface cleaning

### TECHNICAL DESCRIPTION

Product name: Wipes – Embossed Airlaid 60 g/m<sup>2</sup>  
Product type: Single use, non sterile  
Material: Embossed (2 sides) Airlaid nonwoven - cellulosic fibers  
Grammage: 60 g/m<sup>2</sup>  
Origin: France 

### MARKING ON PACKAGING



\* Do not use if package is damaged

Item	Colors	Dimensions (cm)	Unit weight (g)	Packaging
F35000	White	49 x 45	13,2	10 polybags of 50 units
E03302		28 x 30	5,2	30 polybags of 50 units
E03314		28 x 38	6,3	12 polybags of 50 units
E03315	Blue			
E03318	Pink			

## MECHANICAL INFORMATION

TEST	TEST METHOD		AVERAGE	UNIT	TOLERANCE
Grammage	EDANA 40.3-90		57	<i>g/m<sup>2</sup></i>	<i>± 6gsm</i>
Thickness	EDANA 30.5-99		0,9	<i>mm</i>	<i>± 0,1mm</i>
Tensile strength	EDANA 20.2-89	Dry MD :	20	<i>N/5cm</i>	<i>± 2</i>
		Dry CD :	14,4	<i>N/5cm</i>	<i>± 1,5</i>
		Wet MD :	9	<i>N/5cm</i>	<i>± 1</i>
		Wet CD :	7,2	<i>N/5cm</i>	<i>± 1</i>
Absorption	EDANA 10.4-02		10,5	<i>g/g</i>	<i>± 10%</i>

## MICROBIOLOGICAL INFORMATION

Evaluation of the microbial contamination according to standards NF EN ISO 11737 : 2018 + A1 :2021.

## STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.

## LOGISTICAL INFORMATION

OUTER CASE SPECIFICATIONS					
Article	Size (mm)	Weight (kg)	Qty / palet	Qty/ box	EAN code bag
F35000	800 x 400 x 308	7,6	21 boxes (3 boxes x 7 layers)	10 polybags of 50 units	3 662 036 001 993
E03302		9		30 polybags of 50 units	3 662 036 001 313
E03314	600 x 400 x 219	4,5	40 boxes (4 boxes x 10 layers)	12 polybags of 50 units	3 662 036 001 405
E03315					3 662 036 001 412
E03318					3 662 036 001 443