

FICHE TECHNIQUE

WIPES

Spunlace 60 g/m²

28 x 38 cm

Kolmi®



Item	Colour
E09207	<input type="checkbox"/> White

FEATURES & BENEFITS

- Very absorbent and resistant nonwoven sheet

PRECONISATION :

- Used for wiping body
- Used for surface cleaning

TECHNICAL DESCRIPTION

Product name: Wipes – Spunlace 60 g/m²
Product type: Single use, non sterile
Material : Spunlaced nonwoven
50% Viscose fibers / 50% Polyester fibers
Grammage : 60 g/m²
Unit Weight: 6 g (±10%)
Size: 280 x 380 mm (±5mm)
Packaging: Bag of 70 units
Colour : White
Origin : France

MARKING ON PACKAGING



* Do not use if package is damaged

MECHANICAL INFORMATION

TEST	TEST METHOD	AVERAGE	UNIT	TOLERANCE	
Grammage	NWSP 130.1.R0	55	g/m ²	± 10%	
Thickness	NWSP 120.6.R0	0,57	mm	± 0,12mm	
Tensile strength	NWSP 110.4.R0	Dry MD :	85	N/5cm	min 55
		Dry CD :	21	N/5cm	min 14
Elongation		Dry MD :	35	%	max 47
Dry CD :		125	%	max 155	
Absorption	NWSP 010.1.R0	900	%	min 600%	

MICROBIOLOGICAL INFORMATION

Assessment of initial microbial contamination in accordance with NF EN ISO 11737: 2018 + A1: 2021.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTICAL INFORMATION

OUTER CASE SPECIFICATIONS					
Article	Size (mm)	Weight (kg)	Qty/ palet	Qty/ box	EAN code bag
E09207	600 x 400 x 219	5,6	40 boxes (4 boxes x 10 layers)	12 bags of 70 units	3 662 036 001 603