

FICHE TECHNIQUE

WIPES

Plain Airlaid 80 g/m²

120 x 80 cm

Kolmi®



Item	Colour
F35004	<input type="checkbox"/> White


FEATURES & BENEFITS

- Very absorbent and resistant nonwoven sheet

PRECONISATION :

- Used for wiping body

TECHNICAL DESCRIPTION

Product name: Wipes – Plain Airlaid 80 g/m²
Product type: Single use, non sterile
Material : Plain Airlaid nonwoven - cellulosic fibers
Grammage : 80 g/m²
Unit Weight: 77 g (±10%)
Size: 120 x 80 cm (±5mm)
Packaging: Bag of 50 units
Colour : White
Origin : France 

MARKING ON PACKAGING



* Do not use if package is damaged

MECHANICAL INFORMATION

TEST	TEST METHOD		AVERAGE	UNIT	TOLERANCE
Grammage	EDANA 40.3-90		75	g/m ²	± 5%
Thickness	EDANA 30.5-99		1,05	mm	± 0,10mm
Tensile strength	EDANA 20.2.89	Dry MD :	26	N/5cm	± 2,6
		Dry CD :	20	N/5cm	± 2
		Wet MD :	14	N/5cm	± 1,5
		Wet CD :	10	N/5cm	± 1
Absorption	EDANA 10.4-02		11	g/g	± 10%

MICROBIOLOGICAL INFORMATION

Assessment of initial microbial contamination in accordance with NF EN ISO 11737 : 2018 + A1 :2021.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTICAL INFORMATION

OUTER CASE SPECIFICATIONS					
Article	Size (mm)	Weight (kg)	Qty/ palet	Qty/ box	EAN code bag
F35004	800 x 400 x 220	4,5	30 boxes (3 boxes x 10 layers)	1 bag of 50 units	3 662 036 002 020