

# TECHNICAL DATA SHEET

## Safetouch® Vitals™

Wash glove – Needle punch

70g/m<sup>2</sup>

# Kolmi®



### FEATURES AND BENEFITS

- Absorbent, soft and very resistant washing glove
- Ultrasonic sealing
- Compact & hygienic vacuum polybag

### RECOMMENDED USE

- Used for body hygiene

### PRODUCT DESCRIPTION

Product name: Washglove – Needle punch  
Product type: Single use, non-sterile  
Material: Needle punch nonwoven 100% Polyester  
Ultrasonic sealing  
Colour: White  
Packaging: 70 or 50 units in vacuum polybag  
Country of origin: China



Article	Color	Dimensions (cm)	Unit weight (g)	Packaging
G01421	White	14 x 21	4,1	20 bags x 70 units
G01521		15 x 21	4,3	
G01623		16 x 23	5,5	

## PRODUCT PERFORMANCES

Test	Standard		Average	Unit	Tolérance
Grammage	ISO 9073-1		70	<i>g/m<sup>2</sup></i>	$\pm 10\%$
Thickness	ISO 9073-2		0,33	<i>mm</i>	$\pm 10\%$
Tensile Strength	ISO 9073-3	Wet MD:	67	<i>N/5cm</i>	$\pm 10\%$
		Wet CD:	44	<i>N/5cm</i>	$\pm 10\%$
Elongation	ISO 9073-3	Wet MD:	56	%	$\pm 10\%$
		Wet CD:	70	%	$\pm 10\%$
Absorption	ISO 9073-6		500	%	$\pm 20\%$

## STORAGE CONDITIONS

Normal conditions of conservation & storage: must be stored in a cool or room temperature environment and in a dry atmosphere.

## PACKAGING MARKING



\* Do not use if packaging is damaged

## LOGISTICAL INFORMATION

### OUTER CASE SPECIFICATIONS

ITEM	SIZE (mm)	WEIGHT (kg)	QTE/ PALLET	OUTER CASE EAN	QTE / OUTER CASE	POLYBAG EAN
G01421	300 x 220 x 430	6	36 boxes 6 layers of 6	3 662 036 021 243	20 Polybags of 70 units	3 662 036 021 236
G01521	320 x 220 x 320	4,7	60 boxes 5 layers of 12	3 662 036 021 229	20 Polybags of 50 units	3 662 036 021 212
G01623	340 x 240 x 320	5,5	60 boxes 10 layers of 6	3 662 036 019 738		3 662 036 019 745