


FICHE TECHNIQUE

ACTIVE WASH GLOVE

Airlaid + PE 70g/m²

15 x 23 cm

Kolmi®



FEATURES & BENEFITS

- Very absorbent, soft and very résistent washing glove

PRECONISATION :

- Used for body hygiene

TECHNICAL DESCRIPTION

Product name: Washglove – Airlaid/PE
Product type: Single use, non sterile
Material : Embossed Airlaid gluing on blue polyethylene film
Grammage : 70 g/m²
Unit Weight: 4,5 g (±10%)
Size: 150 x 230mm (±5mm)
Packaging: Bag of 50 units
Colour : White / Blue
Origin : China

Item	Colour
G04301	<input type="checkbox"/> White / Blue

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.

MARKING ON PACKAGING



* Do not use if package is damaged

LOGISTICAL INFORMATION

OUTER CASE SPECIFICATIONS						
Article	Size (mm)	Weight (kg)	Qty / palet	Qty/ box	Box EAN code	Bag EAN code
G04301	330 x 240 x 390	4,5	50 boxes (10 boxes x 5 layers)	20 bags of 50 units	3 662 036 011 022	3 662 036 011 039