

FICHE TECHNIQUE

WASHGLOVE

NEEDLE PUNCH 70 g/m²

15 x 22 cm

Kolmi®



FEATURES 1 BENEFITS

- Very absorbent, soft and very résistent washing glove
- Ultrasonic sealing

PRECONISATION :

- Used for body hygiene

TECHNICAL DESCRIPTION

Product name: Washglove – Needle punch
 Product type: Single use, non sterile
 Material : Needle punch nonwoven 100% Polyester
 Unit Weight: 5g (±10%)
 Size: 150 x 220mm (±5mm)
 Packaging: Vacuum polybag of 50 units
 Colour : White
 Origin : China

Item	Colour
G05028	<input type="checkbox"/> White

MARKING ON PACKAGING



* Do not use if package is damaged

PRODUCT PERFORMANCES

Test	Standard		Average	Unit	Tolérance
Grammage	ISO 9073-1		70	g/m ²	± 10%
Thickness	ISO 9073-2		0,33	mm	± 10%
Tensile Strength	ISO 9073-3	Dry MD:	96	N/5cm	± 10%
		Dry CD:	78	N/5cm	± 10%
		Wet MD:	114	N/5cm	± 10%
		Wet CD:	68	N/5cm	± 10%
Elongation	ISO 9073-3	Dry MD:	54	%	± 10%
Absorption	ISO 9073-6		600	%	± 20%

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTICAL INFORMATION

OUTER CASE SPECIFICATIONS						INNER CASE SPECIFICATIONS
Article	Size (mm)	Weight (kg)	Qty / palet	Qty/ box	OUTERCASE EAN	POLYBAG EAN
G05028	320 X 230 X 320	4,5	60 boxes (6 boxes x 10 layers)	20 vacuum bags of 50 units	3 662 036 011 947	3 662 036 021 250