

# FICHE TECHNIQUE

## ACTIVE WASH GLOVE

Spunlace honeycomb 80g/m<sup>2</sup>

15 x 22 cm

# Kolmi®



Item	Colour
G09010	<input type="checkbox"/> White

### FEATURES 1 BENEFITS

- Very absorbent, soft and very résistent washing glove
- Ultrasonic sealing

### PRECONISATION :

- Used for body hygiene

### TECHNICAL DESCRIPTION

Product name:	Washglove – Honeycomb Spunlace
Product type:	Single use, non sterile
Material :	Spunlace
Grammage :	80 g/m <sup>2</sup>
Unit Weight:	5g (±10%)
Size:	150 x 220mm (±5mm)
Packaging:	Bag of 50 units
Colour :	White
Origin :	China

### MECHANICAL INFORMATION

TEST	STANDARD		AVERAGE	UNIT	TOLERANCE
Grammage	ISO 9073-1		80	<i>g/m<sup>2</sup></i>	<i>± 5%</i>
Thickness	ISO 9073-2		0,5	<i>mm</i>	<i>± 10%</i>
Tensile strength	ISO 9073-3	Dry MD :	169	<i>N/5cm</i>	<i>± 10%</i>
		Dry CD :	119	<i>N/5cm</i>	<i>± 10%</i>
		Wet MD :	136	<i>N/5cm</i>	<i>± 10%</i>
		Wet CD :	124	<i>N/5cm</i>	<i>± 10%</i>
Elongation	ISO 9073-3		60	<i>%</i>	<i>± 10%</i>
Absorption	ISO 9073-6		800	<i>%</i>	<i>± 10%</i>

## MICROBIOLOGICAL INFORMATION

Assessment of initial microbial contamination in accordance with NF EN ISO 11737 : 2018 + A1 :2021.

## STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.

## MARKING ON PACKAGING



\* Do not use if package is damaged

## LOGISTICAL INFORMATION

OUTER CASE SPECIFICATIONS						
Article	Size (mm)	Weight (kg)	Qty / palet	EAN code box	Qty/ box	EAN code bag
G09010	470 x 310 x 360	5,5	25 boxes (5 boxes x 5 layers)	3 700 381 511 498	20 bags of 50 units	3 700 381 511 504