

# Technical Data Sheet

## SAFE FEET® SKIDguard

Universal Shoe Cover  
Water proof White Sole  
Medical device Class I



### FEATURES & BENEFITS

- Nonwoven spunbonded 100% polypropylene.
- Gluing sealing.
- Elastic free of latex at the ankle.
- Supple and very resistant nonwoven.
- Water proof and anti-slip white sole 60µm.
- They adapt to different morphologies and all types of shoes.

### RECOMMENDED USE

- Can be used for staff or patients in operating theatres and care services.
- Used in industrial environment and food industry.

### PRODUCT DETAILS

Product name : Shoe Cover White PP SafeFeet® SkidGuard™  
 Product type : Single use, non sterile  
 Material : Nonwoven Polypropylene Spunbond + Polyethylene Sole  
 Confection : Gluing  
 Color : White + White sole  
 Size : 420 mm x 230 mm (± 10%)  
 Modele XXL + Very extra high  
 Unit Weight: 14,4g (± 10%)  
 Packaging : Dispenser box of 300 pcs  
 Manufacturing country : FRANCE ■ ■

Item	Color
H07101A	<input type="checkbox"/> White

### MECHANICAL INFORMATION

Measurements done by independent laboratories  
(measures are available upon request)

TEST	STANDARD	LABORATORY		AVERAGE	UNIT	TOLERANCE
Grammage	ISO 9073-1	Cetelor		35,00	g/m <sup>2</sup>	±10%
Tensile strength	ISO 9073-3	Cetelor	Tensile stren. MD:	64,90	N	±10%
			Tensile stren. CD:	80,20	N	±10%
Elongation			Elongation MD:	50,74	%	±10%
			Elongation CD:	63,24	%	±10%

## MICROBIOLOGICAL INFORMATION

Assessment of initial microbial contamination according to standards NF EN ISO 11737:2018 + A1: 2021.  
Additional microbiological controls: ASR, E.coli, Staphylococci.

## CERTIFICATIONS & STANDARDS

Meet the requirements of Regulation (EU) 2017/745 on medical devices.  
Complies with applicable standards EN ISO 14971.  
Manufacturing site : ISO 9001 – ISO 13485.

## PRECAUTIONS

The product must be used on intact skin only.  
The reuse or extended use may lead to infection or cross-contamination.  
Please refer to your national legislation on the subject of disposal of this product once it has been used.

## STORAGE

Normal conditions of conservation & storage : not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.  
Shelf life of the product: 5 years.

## MARKING ON PACKAGING



\* Do not use if package is damaged.

## LOGISTICS INFORMATION

OUTER CASE SPECIFICATIONS				
ITEM	SIZE (mm)	GROSS WEIGHT (kg)	QTY / PALLET	BARCODE
H07101A	400 x 400 x 370	4,9	36 cases	3 662 036 002 730