

TECHNICAL DATA SHEET

ISOLATION GOWN POLYPROPYLENE



FEATURES & BENEFITS

- Comfortable and light, the material is highly breathable
- Without collar and elasticated wrists
- Back closure with laces at the belt and the neck
- Without pocket
- Long sleeves
- Single size fits all

Recommended use:

- Worn by caregivers and visitors over clothing, it protects patients from the risks of external infection



PRODUCT DETAILS

Product name :	Isolation gowns
Product type :	Single use, non sterile
Material :	Polypropylene
Confection :	Sewn, without collar, without pocket, long sleeves, elasticated wrists, Back closure with laces at the belt and the neck.
Grammage :	17 g/m ² (+/-3)
Length :	105 cm (+/- 1cm)
Width :	134 cm (+/- 1cm)
Unit weight :	30 g (+/- 10%)
Packaging :	10 polybags of 10 units
Manufacturing country:	China

Item	Color
J40451	Green
J40452	White
J40453	Sky blue
J40454	Yellow
J40456	Pink

CERTIFICATIONS & STANDARDS

Manufactured under ISO 9001 certified quality system

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTICS INFORMATION



Outer case

Item	Size Mm	Gross Weight kg	QTY/ pallet	QTY / case
J40451 J40452 J40453 J40454	460 x 280 x 260	3,7	36 cases (6 layers of 6 cartons)	10 polybags of 10 pieces

EAN CODES

Item	Outer case	Polybag
J40451	3 662 036 019 271	3 662 036 019 288
J40452	3 662 036 019 295	3 662 036 019 301
J40453	3 662 036 019 318	3 662 036 019 325
J40454	3 662 036 019 332	3 662 036 019 349