

TECHNICAL DATA SHEET

PATIENT GOWN POLYPROPYLENE



FEATURES & BENEFITS

- Thanks to its higher grammage than isolation gowns, patient gowns permit better opacity
- Pleasant and light, the polypropylene is very breathable
- Collarless model
- Lace up back closure
- One size fits all
- Individually packaged for hygiene and controlled distribution
- No pockets to prevent objects from falling out

PRODUCT DETAILS

Product name:	Patient Gown
Product type:	Single use, non-sterile
Material:	Polypropylene
Confection:	Sewn, no collar, no pocket
Colour:	Navy blue
Packaging:	Polybag of 10 units
Origin:	China



J40459

J40460

Item	Sleeves	Back closure with laces	Grammage (±10%)	Length (±1 cm)	Width (±1 cm)	Unit weight (±10%)
J40459	Short sleeves	At the belt	34 g/m ²	110 cm	140 cm	53 g
J40460	3/4 sleeves	At the belt and the neck	42 g/m ²	115 cm	140 cm	84 g

CERTIFICATIONS & STANDARDS

Manufactured under ISO 9001 certified quality system

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTICS INFORMATION



Outer case

Item	Size Mm	Gross Weight Kg	QTY/ pallet	QTY / case
J40459	460 x 280 x 260	6,1	30 cases (5 layers of 6 cases)	10 polybags of 10 pieces
J40460	480 x 290 x 420	9,3	30 cases (5 layers of 6 cases)	

EAN CODES

Item	Outer case	Polybag
J40459	3 662 036 019 370	3 662 036 019 387
J40460	3 662 036 019 714	3 662 036 019 721