

TECHNICAL DATA SHEET

PP CLEAN AIR SUIT FOR HEALTHCARE WORKERS



FEATURES & BENEFITS

- Comfortable and light, the material is highly breathable
- The clean air suit provides optimum coverage of the upper body.
- The ample shape allows the user to have a good comfort
- The white jersey cuffs and collar guarantee better comfort.
- Individually wrapped for hygiene and controlled distribution.

PRODUCT DETAILS

Product name :	Clean air suit for Healthcare workers
Product type :	Single use, non sterile
Material :	Polypropylene
Confection :	Sewn, jersey collar, elasticated jersey cuffs, two pockets, snap fastener
Color :	Dark blue
Grammage :	40 g/m ² (+/-10%)
Length :	78 cm (+/-1cm)
Width :	68 cm (+/-1cm)
Packaging :	50 individual polybags
Manufacturing country :	China

CERTIFICATIONS & STANDARDS

Manufactured under ISO 9001 certified quality system

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

INFORMATIONS LOGISTIQUES

Outercases

Article	Taille mm	Poids brut kg	QTE/ palette	CODE EAN
J40895/L	460 x 300 x 300	6,0	36 outercases (6 layers of 6 outercases)	3 662 036 021 434

Polybag

QTE / carton	CODE EAN
50 individual polybags	3 662 036 021 441