

TECHNICAL DATA SHEET

CPE APRON GOWN SPLASH RESISTANT J50656



FEATURES & BENEFITS

- Provides superior protection to a traditional apron
- Covers almost the whole body but also the arms
- Resistant PE tie bonds.
- « Thumb loops » system to avoid the sleeves go back
- It adapts to different morphologies

PRODUCT DETAILS

Product name :	CPE apron gown splash resistant
Product type :	Single use
Material :	CPE
Confection :	Cutting, folding, sewing
Color :	Light blue
Thickness :	35 µm
Dimensions :	119 x 105 cm
Unit weight :	54 g
Manufacturing country :	China

CERTIFICATIONS & STANDARDS

Manufactured under ISO 9001 qualified system.

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTICS INFORMATION

OUTERCASE

Item	Size mm	Gross weight kg	QTY/pallet	CODE EAN
J50656	380 x 200 x 310	3,9	48 cases	3 662 036 019 961

POLYBAG

QTY / case	CODE EAN
3 bags of 25 pcs	3 662 036 019 707