

## MBKN3 NITRILE GLOVE SAFETOUCH® ADVANCED™ BLACK LITE

Class I Medical Device and Category III Personal Protective Equipment, Powder Free

### FEATURES AND BENEFITS

- Great elasticity.
- Good dexterity.
- Its smooth appearance gives it correct touch and grip of objects.
- Ultra comfortable, this generation allows an excellent stretchability and an optimal comfort.
- Watertight and airtight.

### Recommendations:

- Glove for professional use.
- In the healthcare environment, this glove is used to make barrier against microorganisms between the professional and the patient and their environment.
- In an industrial environment, this glove is intended to protect the wearer from tested chemicals and microorganisms (bacteria, moulds and viruses).
- The powder-free glove reduces the risk of airborne contamination and limits the risk of irritation.

### PRODUCT DETAILS

Product name:	<b>SAFETOUCH ADVANCED Black Lite</b>
Product type:	Single use, no sterile
Material:	Nitrile
Colour:	Black
Shape:	Ambidextrous
Cuff:	Rolled
Surface:	Textured fingertips
Powdered:	No
Palm thickness:	0,07 (±0,02mm)
Fingers thickness:	0,10 (± 0,02mm)
Glove length:	≥ 240 mm
Unit weight (M):	3,5g (± 0,2)
Manufacturing country:	Malaysia

Item	Size
MBKN3-6XS	T 6 XS
MBKN3-7S	T 7 S
MBKN3-8M	T 8 M
MBKN3-9L	T 9 L
MBKN3-10XL	T 10 XL

Marking on packaging:



### INDUSTRIES

Health / Medical / Dental / Industry / Food

## CERTIFICATION AND STANDARDS

Manufactured under an ISO 13485 and ISO 9001 certified quality system.

In accordance with requirements of Reg (EU) 2017/745, UK Medical Device Regulation 2002 and the Swiss Ordonnance (ODim) from July 1st 2020, on medical devices.

Comply with the standards EN 455-1/2/3/4 for single use medical gloves.

In accordance with requirements of Reg (EU) 2016/425 and Regulation 2016/425 as amended to apply in GB, on personal protective equipments.

Comply with the EN 420, EN ISO 374-1/2/4/5, for protective against tested chemicals and microorganisms including viruses.

In accordance with requirements of Reg.(EC) 1935/2004 and the French decree of 5th August, 2020. This glove is suitable for food contact, with all types of foodstuff (see annex 1).

## PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.

The reuse or extended use may lead to infection or cross-contamination.

Please refer to your national legislation on the subject of disposal of this product once it has been used.

## STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, store in a cool environment or at room temperature.

Shelf life of the product: 3 years.

## LOGISTIC INFORMATION








### Specifications outer case



### Specifications inner case

### Color code inner and outer

Item	Size mm	Gross Weight kg	QTY/pallet	Size mm	QTY	Color
MBKN3-6XS		4,1				
MBKN3-7S		4,7				
MBKN3-8M	262 x 252 x 257	5,1	60 cases (6 layers of 10 cases)	245 x 125 x 48	10 dispenser boxes of 100 units	
MBKN3-9L		5,3				
MBKN3-10XL		6,2				



### EAN Codes

Item	Outer case	Inner case
MBKN3-6XS	3 662 036 020 802	3 662 036 020 857
MBKN3-7S	3 662 036 020 819	3 662 036 020 864
MBKN3-8M	3 662 036 020 826	3 662 036 020 871
MBKN3-9L	3 662 036 020 833	3 662 036 020 888
MBKN3-10XL	3 662 036 020 840	3 662 036 020 895