

MCLV5
SAFETOUCH EVERSTRONG VINYL POWDER-FREE GLOVE
 Category III Personal Protective Equipment & Class I Medical Device

FEATURES & BENEFITS

- Suitable elasticity
- Good dexterity
- Good feeling when touching objects and good grip on objects by its smooth appearance
- Watertight and airtight

Recommendation:

- Glove for professional use
- In the healthcare environment, this glove is used to make barrier against microorganisms between the healthcare professional, the patient and their environment
- In an industrial environment, this glove is used to protect the wearer from microorganisms (bacteria, molds and viruses) and tested chemicals
- For food handling / preparation (HACCP)

PRODUCT DESCRIPTION

Product name:	SafeTouch Everstrong
Product type :	Single use, non-sterile
Material :	Vinyl
Colour :	Clear
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Smooth
Powdered :	No
Palm thickness :	0.08 (±0,02 mm)
Finger thickness :	0.08 (±0,02 mm)
Glove length:	240mm minimum
Unit weight (M):	5,3 (±0,2 g)
Country of origin :	China

Item	Size	
MCLV5-7S	T 7	S
MCLV5-8M	T 8	M
MCLV5-9L	T 9	L
MCLV5-10XL	T 10	XL

MARKING ON PACKAGING

INDUSTRIES

Food / Non residential catering / Community / Administration / Tatoon artist/ Sanitary / Maintenance / Cleaning / Laboratory / Pharmaceutical / Research

CERTIFICATION

Manufactured under an ISO 13485 and ISO 9001 certified quality system

In accordance with requirements of Reg (UE) 2017/745 for medical device and UK Medical Device Regulation 2002.

Comply with the EN 455-1/2/3/4 for single use medical gloves for single use.

In accordance with requirements of Reg (UE) 2016/425 for personal protective equipment and Regulation 2016/425 on personal protective equipment, as amended to apply in GB.

Comply with the EN ISO 21420, EN ISO 374-1/2/4/5 for protection against tested chemicals and microorganisms including viruses.

In accordance with requirements of Reg.(CE) 1935/2004 and Reg.(EU) 10/2011. This glove is suitable for food contact with all food*

* Restriction : It is not recommended to put vinyl powdered gloves in contact with foods containing fats or oils (see Annex 1).

PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.

The reuse or extended use may lead to infection or cross-contamination.

Please refer to your national legislation on the subject of disposal of this product once it has been used.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at room temperature.

Shelf life of the product: 3 years

LOGISTIC INFORMATION







Specifications outer case



Specifications inner case

Colour code inner and outer

Item	Size mm	Gross weight kg	QTY/ Pallet	Size mm	QTY	Colour
MCLV5-7S	315 x 259 x 245	5,18	60	230 x 123 x 60	10 x 100	
MCLV5-8M	315 x 259 x 245	6,03	60	230 x 123 x 60	10 x 100	
MCLV5-9L	315 x 259 x 245	6,36	60	230 x 123 x 60	10 x 100	
MCLV5-10XL	315 x 259 x 245	6,90	60	230 x 123 x 60	10 x 100	



EAN Codes

ITEM	Inner case	Outer case
MCLV5-7S	3 700 381 522 012	3 700 381 522 005
MCLV5-8M	3 700 381 522 036	3 700 381 522 029
MCLV5-9L	3 700 381 522 050	3 700 381 522 043
MCLV5-10XL	3 700 381 522 074	3 700 381 522 067