

**MCLVP4**
**SAFETOUCH EVERSTRONG VINYL POWDERED GLOVE**  
 Category III Personal Protective Equipment

**FEATURES & BENEFITS**

- Suitable elasticity
- Good dexterity
- Good feeling when touching objects and good grip on objects by its smooth appearance
- Watertight and airtight


**Recommendation:**

- Glove for professional use
- This glove is used to protect the wearer from microorganisms (bacteria, molds and viruses) and tested chemicals
- For food handling / preparation (HACCP)

**PRODUCT DESCRIPTION**

Product name:	<b>SafeTouch EverStrong</b>
Product type :	Single use, non-sterile
Material :	Vinyl
Colour :	Clear
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Smooth
Powdered :	Yes
Palm thickness :	0.08 (±0,02 mm)
Finger thickness :	0.08 (±0,02 mm)
Glove length:	240mm minimum
Unit weight (M):	4,5 (±0,2 g)
Country of origin :	China

Item	Size	
MCLVP4-7S	T 7	S
MCLVP4-8M	T 8	M
MCLVP4-9L	T 9	L
MCLVP4-10XL	T 10	XL

**MARKING ON PACKAGING**

**INDUSTRIES**

Food / Non residential catering / Community / Administration / Sanitary / Maintenance / Cleaning / Laboratory / Pharmaceutical / Research

### CERTIFICATION

Manufactured under an ISO 13485 and ISO 9001 certified quality system.

In accordance with requirements of Reg (UE) 2016/425 for personal protective equipment and Regulation 2016/425 on personal protective equipment, as amended to apply in GB.

Comply with the EN ISO 21420, EN ISO 374-1/2/4/5 for protection against tested chemicals and microorganisms including viruses.

In accordance with requirements of Reg.(CE) 1935/2004 and Reg.(EU) 10/2011. This glove is suitable for food contact with all food\*

\* Restriction : It is not recommended to put vinyl powdered gloves in contact with foods containing fats or oils (see Annex 1).

### PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.

The reuse or extended use may lead to infection or cross-contamination.

Please refer to your national legislation on the subject of disposal of this product once it has been used.

### STORAGE CONDITIONS

Normal conditions of conversation & storage: not to be exposed to moisture and sun, must be stored at room temperature.

Shelf life of the product: 3 years

### LOGISTIC INFORMATION







#### Specifications outer case



#### Specifications inner case

#### Colour code inner and outer

Item	Size mm	Gross weight kg	QTY/ Pallet	Size mm	QTY	Colour
MCLVP4-7S	315 x 259 x 245	5,15	60	230 x 123 x 60	10 x 100	
MCLVP4-8M	315 x 259 x 245	5,87	60	230 x 123 x 60	10 x 100	
MCLVP4-9L	315 x 259 x 245	6,63	60	230 x 123 x 60	10 x 100	
MCLVP4-10XL	315 x 259 x 245	7,02	60	230 x 123 x 60	10 x 100	



#### EAN Codes

ITEM	Inner case	Outer case
MCLVP4-7S	3 700 381 522 098	3 700 381 522 081
MCLVP4-8M	3 700 381 522 111	3 700 381 522 104
MCLVP4-9L	3 700 381 522 135	3 700 381 522 128
MCLVP4-10XL	3 700 381 522 159	3 700 381 522 142