

1136 SAFETOUCH 400 POWDER-FREE NITRILE GLOVE

Class I Medical Device and Category III Personal Protective Equipment

FEATURES & BENEFITS

- 400 mm length for extra protection.
- Polymer coated for easy donning.
- Its textured appearance gives it correct touch and grip of objects.
- Ultra comfortable, this glove allows an excellent stretchability and an optimal comfort.

Recommendations:

- Glove for professional use.
- In the healthcare environment, this glove acts as a barrier against microorganisms between the professional and the patient and their environment.
- In an industrial environment, this glove is intended to protect the wearer from tested chemicals (including chemotherapy medicines) and microorganisms (bacteria, moulds and viruses).
- The powder-free glove reduces the risk of airborne contamination and limits the risk of irritation.
- For food handling / preparation (HACCP).




PRODUCT DESCRIPTION

Product name :	SafeTouch 400
Product type :	Single use, non-sterile
Material :	Nitrile
Colour :	Blue
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Textured
Powdered :	No
Palm thickness :	0,15 (±0,02mm)
Finger thickness :	0,16 (±0,02mm)
Glove Length:	400mm
Unit weight (M):	14,5g (± 0,5g)

Article	Size	
1136-A	T6	XS*
1136-B	T7	S
1136-C	T8	M
1136-D	T9	L
1136-E	T10	XL

Marking on packaging:



INDUSTRIES

Medical/ Food/ Laboratory/ Industry

CERTIFICATION

Manufactured under an ISO 13485 and ISO 9001 certified quality system.
 In accordance with requirements of Reg (EU) 2017/745, on medical devices.
 Comply with the standards EN 455-1/2/3/4 for single use medical gloves.
 In accordance with requirements of Reg (EU) 2016/425, on personal protective equipments.
 Comply with the EN ISO 21420, EN ISO 374-1/2/4/5 and ASTM D6978-05, for protective gloves against tested chemicals and microorganisms including viruses.
 In accordance with requirements of Reg.(EC) 1935/2004 and the French decree of 5th August 2020. This glove is suitable for food contact, with all types of foodstuff.

PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.
 Single use only, do not reuse.
 Please refer to your national legislation about disposal of this product once it has been used.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, store in a cool environment or at room temperature.
 Shelf life of the product: 3 years.

LOGISTIC INFORMATION

**Specifications
outer case**



**Specifications
inner case**



**Color code
inner and outer**

Item	Size mm	Gross weight Kg	QTY/ pallet	Size mm	QTY	Color
1136-A		7				
1136-B		7,75				
1136-C	415 x 275 x 250	8,75	36 cartons (6 layers of 6 cartons)	400 x 130 x 45	10 boxes of 50 pieces	
1136-D		9,25				
1136-E		10				

EAN Codes

Item	Outer case	Inner case
1136-A	3701407401052	3701407401069
1136-B	3701407400567	3701407400550
1136-C	3701407400581	3701407400574
1136-D	3701407400604	3701407400598
1136-E	3701407400628	3701407400611

