

FICHE TECHNIQUE

Kolmi®

WIPES

Plain Airlaid 80 g/m²

118 x 75 cm



FEATURES & BENEFITS


- Very absorbent and resistant nonwoven sheet

PRECONISATION :

- Used for wiping body

Item	Colour
F35004	<input type="checkbox"/> White

TECHNICAL DESCRIPTION

Product name:	Wipes – Plain Airlaid 80 g/m ²
Product type:	Single use, non sterile
Material :	Plain Airlaid nonwoven - cellulosic fibers
Grammage :	80 g/m ²
Unit Weight:	70 g (±10%)
Size:	118 x 75 cm (±2,5%)
Packaging:	Bag of 50 units
Colour :	White
Origin :	 Europe

MARKING ON PACKAGING



* Do not use if package is damaged

MECHANICAL INFORMATION

TEST	TEST METHOD		AVERAGE	UNIT	TOLERANCE
Grammage	EDANA 40.3-90		75	g/m ²	± 5%
Thickness	EDANA 30.5-99		1,05	mm	± 0,10mm
Tensile strength	EDANA 20.2.89	Dry MD :	26	N/5cm	± 2,6
		Dry CD :	20	N/5cm	± 2
		Wet MD :	14	N/5cm	± 1,5
		Wet CD :	10	N/5cm	± 1
Absorption	EDANA 10.4-02		11	g/g	± 10%

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTICAL INFORMATION

OUTER CASE SPECIFICATIONS					
Article	Size (mm)	Weight (kg)	Qty/ palet	Qty/ box	EAN code bag
F35004	590 x 410 x 215	4,3	32 boxes (4 boxes x 8 layers)	1 bag of 50 units	3 662 036 002 020