

MCLV5

SAFETOUCH EVERSTRONG VINYL POWDER-FREE GLOVE Category III Personal Protective Equipment & Class I Medical Device

FEATURES & BENEFITS

- Suitable elasticity.
- Good dexterity.
- Good feeling when touching objects and good grip on objects by its smooth appearance.
- Watertight and airtight.

Recommendation:

- Glove for professional use.
- In the healthcare environment, this glove is used to make barrier against microorganisms between the healthcare professional, the patient and their environment.
- In an industrial environment, this glove is used to protect the wearer from microorganisms (bacteria, molds and viruses) and tested chemicals .
- For food handling / preparation (HACCP).

PRODUCT DESCRIPTION

| | |
|---------------------|-----------------------------|
| Product name: | SafeTouch Everstrong |
| Product type : | Single use, non-sterile |
| Material : | Vinyl |
| Colour : | Clear |
| Shape : | Ambidextrous |
| Cuff : | Rolled |
| Surface : | Smooth |
| Powdered : | No |
| Palm thickness : | 0.08 (±0,02 mm) |
| Finger thickness : | 0.08 (±0,02 mm) |
| Glove length: | 240mm minimum |
| Unit weight (M): | 5,3 (±0,2 g) |
| Country of origin : | China |

| Item | Size |
|------------|---------|
| MCLV5-7S | T 7 S |
| MCLV5-8M | T 8 M |
| MCLV5-9L | T 9 L |
| MCLV5-10XL | T 10 XL |

MARKING ON PACKAGING



INDUSTRIES

Food / Non residential catering / Community / Administration / Tatoo artist/ Sanitary / Maintenance / Cleaning / Laboratory / Pharmaceutical / Research

CERTIFICATION

Manufactured under an ISO 13485 and ISO 9001 certified quality system.
 In accordance with requirements of Reg (EU) 2017/745, on medical devices.
 Comply with the standards EN 455-1/2/3/4 for single use medical gloves.
 In accordance with requirements of Reg (EU) 2016/425, on personal protective equipments.
 Comply with the EN ISO 21420, EN ISO 374-1/2/4/5, for protective gloves against tested chemicals and microorganisms including viruses.
 In accordance with requirements of Reg.(CE) 1935/2004 and Reg.(EU) 10/2011. This glove is suitable for food contact with all food.*

* Restriction : It is not recommended to put vinyl powdered gloves in contact with foods containing fats or oils (see Annex 1).

PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.
 Single use, do not reuse.
 Please refer to your national legislation about disposal of this product once it has been used.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, store in a cool environment or at room temperature.
 Shelf life of the product: 3 years.

LOGISTIC INFORMATION







Specifications outer case



Specifications inner case

Colour code inner and outer

| Item | Size mm | Gross weight kg | QTY/ Pallet | Size mm | QTY | Colour |
|------------|-----------------|-----------------|----------------|----------------|----------|---|
| MCLV5-7S | 315 x 250 x 230 | 5,18 | 60 | 220 x 120 x 60 | 10 x 100 |  |
| MCLV5-8M | 315 x 250 x 230 | 6,03 | 60 | 220 x 120 x 60 | 10 x 100 |  |
| MCLV5-9L | 315 x 250 x 230 | 6,36 | 60 | 220 x 120 x 60 | 10 x 100 |  |
| MCLV5-10XL | 315 x 250 x 230 | 6,90 | 60 | 220 x 120 x 60 | 10 x 100 |  |

EAN Codes

| ITEM | Inner case | Outer case |
|------------|-------------------|-------------------|
| MCLV5-7S | 3 700 381 522 012 | 3 700 381 522 005 |
| MCLV5-8M | 3 700 381 522 036 | 3 700 381 522 029 |
| MCLV5-9L | 3 700 381 522 050 | 3 700 381 522 043 |
| MCLV5-10XL | 3 700 381 522 074 | 3 700 381 522 067 |