


# Paper Mouth Rinse Cups

## FEATURES & BENEFITS

- Paper cups are used to issue a measured amount of liquid (water or oral disinfectant) for rinsing patient's mouth
- Suitable for most common dispensers
- Good resistance to heat and cold thanks to a grammage of 190g/m<sup>2</sup>

## PRODUCT DESCRIPTION

Product name: Paper cups  
Product type: Single use, non-sterile  
Material: Paper  
Grammage: 190g/m<sup>2</sup>  
Capacity: 180ml  
Size: Height: 8,0 cm (± 5%)  
Top diameter: 7,0 cm (± 5%)  
Bottom diameter: 4,5 cm (± 5%)  
Unit weight: 2,9 g (± 10%)  
Packaging: 3.000 pcs  
Colours:



Country of origin: China  
For dental use only.

## INDUSTRIES

Medical / Dental



## Product

Medical Device, Class I / Additional information

## TECHNICAL DATA SHEET

English

### STORAGE CONDITIONS

Normal conditions of conservation & storage:  
not to be exposed to moisture and sun.

### LOGISTICAL INFORMATION

#### Outer case specifications

REF	Colour	Size (mm)	Gross weight (kg)	Content Qty/outer case	EAN Code
120	Blue	720 x 440 x 360	10,0	60 polybags of 50 pcs.	3700381525785
121	White	720 x 440 x 360	10,0	60 polybags of 50 pcs.	3700381525792
122	Green	720 x 440 x 360	10,0	60 polybags of 50 pcs.	3700381525808
123	Yellow	720 x 440 x 360	10,0	60 polybags of 50 pcs.	3700381525815
124	Lavender	720 x 440 x 360	10,0	60 polybags of 50 pcs.	3700381525822
125	Pink	720 x 440 x 360	10,0	60 polybags of 50 pcs.	3700381525839
126	Orange	720 x 440 x 360	10,0	60 polybags of 50 pcs.	3700381525846
127	Dark Blue	720 x 440 x 360	10,0	60 polybags of 50 pcs.	3700381525853
128	Color Mix	720 x 440 x 360	10,0	10 polybags of 50 pcs. per each of the colours blue, green, lavender, pink, orange, dark blue (total 60 polybags of 50 pieces)	3700381525860

