

30.505 SafeBasics™ PE Glove – Translucent orange – 920mm

FEATURES & BENEFITS

- Its polyethylene raw material offers a good dexterity.
- Ambidextrous powered free glove.

TECHNICAL DESCRIPTION

Product name:	SafeBasics™ PE Glove
Product type:	Single use, non sterile
Material :	Polyethylene
Color :	Translucent orange
Shape :	Ambidextrous
Powdered :	Powder free
Texture :	Embossed
Unit weight (M):	9 g
Thickness :	25μ
Origin :	Asia
Packaging:	20 polybags of 100 units



STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISITCAL INFORMATION



Specifications Outer case

Item	Size Mm	Gross Weight (kg)	QTY	QTY/pallet	EAN Code
30.505	390 x 250 x 260	16	20 x 100pcs	20	3 700 381 507 620

Specifications Polybag

Size Mm	Gross Weight (kg)
330 x 230	0,828