

10.110

100 % PE SHOECOVERS

Thickness : 80 µm

FEATURES & BENEFITS

- Extremely durable and fully waterproof shoe covers, thanks to their 100% polyethylene material and 80 µm thickness.
- Suitable for outdoor use.
- Particularly appreciated for use on wet or rough surfaces.
- The material provides good mechanical resistance to abrasion as well as to certain chemicals, although this model is not considered PPE.

PRODUCT DETAILS

Product name :	100% PE SafeWear shoe covers
Product type :	Single use, non sterile
Composition :	100 % Polyethylene
Confection :	Welded
Thickness :	80 µm
Size :	190 x 380 mm
Unit weight :	13 g
Packaging :	Case with bag of 400 pcs
Origin :	China



CERTIFICATIONS & STANDARDS

Manufacturing site : ISO 9001

STORAGE CONDITIONS

Normal storage and preservation conditions: must not be exposed to moisture or sunlight and should be stored in a cool environment or at room temperature.

LOGISTICS

 Outer case spec.

Item	Packaging type	Qty pcs	Size (mm)	Gross weight (kg)	QTY/pallet
10.110	Carton avec sache	400	400 x 360 x 290	5,4	23

 EAN code

Item	Outer case
10.110	3 700 381 500 249