

1191N NITRILE GLOVE SAFETOUCH VITALS POWDER-FREE

Class I Medical Device and Category III Personal Protective Equipment

FEATURES & BENEFITS

- Latex-like feel.
- Agreeable dexterity.
- Excellent grip in both wet and dry conditions.
- Watertight and airtight.
- Convenient packaging that holds **200 gloves** per box reducing storage space, transportation costs and packaging waste.



RECOMMENDATION:

- Users allergic to latex.
- Glove for professional use.
- In the healthcare environment, this glove is used to make barrier against microorganisms between the healthcare professional, the patient and their environment.
- In an industrial environment, this glove is designed to protect the wearer from tested chemicals and microorganisms (bacteria, molds and viruses).
- The powder-free glove reduces the risk of airborne contamination and limits the risk of irritation.
- For food handling / preparations (HACCP).

PRODUCT DESCRIPTION

| | |
|--------------------|--------------------------|
| Product name: | SafeTouch® Vitals |
| Product type : | Single use, non-sterile |
| Material : | Nitrile |
| Colour : | Blue |
| Shape : | Ambidextrous |
| Cuff: | Rolled |
| Surface : | Textured fingertips |
| Powdered : | No |
| Length : | 240 mm |
| Palm thickness : | 0,07 mm (±0,02mm) |
| Finger thickness : | 0,09 mm (±0,02mm) |
| Unit weight (M): | 3,2 g |

| Item | Size |
|---------|---------|
| 1191N-A | T 6 XS |
| 1191N-B | T 7 S |
| 1191N-C | T 8 M |
| 1191N-D | T 9 L |
| 1191N-E | T 10 XL |

Marking on packaging:



INDUSTRIES

Food / Health / Medical / Dental/ Administration / Hygiene/ Industry

CERTIFICATION

Manufactured under an ISO 13485 and ISO 9001 certified quality system.

In accordance with requirements of Reg (EU) 2017/745, on medical devices.

Comply with the standards EN 455-1/2/3/4 for single use medical gloves.

In accordance with requirements of Reg (EU) 2016/425, on personal protective equipments.

Comply with the EN ISO 21420, EN ISO 374-1/2/4/5, for protective gloves against tested chemicals and microorganisms including viruses.

In accordance with requirements of Reg.(EC) 1935/2004 and the French decree of 5th August 2020. This glove is suitable for food contact, with all types of foodstuff.

PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.

The reuse or extended use may lead to infection or cross-contamination.

Please refer to your national legislation about disposal of this product once it has been used.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored in a cool environment or at room temperature.

Shelf life of the product: 5 years.

LOGISTIC INFORMATION








Specifications outer case



Specifications inner case

Colour code inner and outer

| Item | Size mm | Gross weight kg | QTY/ pallet | Size mm | QTY | Colour |
|---------|-----------------|-----------------|----------------|----------------|----------|---|
| 1191N-A | 415 x 261 x 260 | 6,2 | 42 | 243 x 125 x 81 | 10 x 200 |  |
| 1191N-B | 415 x 261 x 260 | 6,5 | 42 | 243 x 125 x 81 | 10 x 200 |  |
| 1191N-C | 415 x 261 x 260 | 7,1 | 42 | 243 x 125 x 81 | 10 x 200 |  |
| 1191N-D | 415 x 261 x 260 | 8,4 | 42 | 243 x 125 x 81 | 10 x 200 |  |
| 1191N-E | 415 x 261 x 260 | 10,3 | 42 | 243 x 125 x 81 | 10 x 200 |  |

EAN Codes

| Item | Inner case | Outer case |
|---------|-------------------|-------------------|
| 1191N-A | 3 700 381 519 890 | 3 700 381 519 883 |
| 1191N-B | 3 700 381 519 913 | 3 700 381 519 906 |
| 1191N-C | 3 700 381 519 937 | 3 700 381 519 920 |
| 1191N-D | 3 700 381 519 951 | 3 700 381 519 944 |
| 1191N-E | 3 700 381 519 975 | 3 700 381 519 968 |