

## 1193N

### NITRILE GLOVE SAFETOUCH VITALS POWDER-FREE

Class I Medical Device and Category III Personal Protective Equipment

#### FEATURES & BENEFITS

- Latex-like feel.
- Agreeable dexterity.
- Excellent grip in both wet and dry conditions.
- Watertight and airtight.



#### RECOMMENDATION:

- Users allergic to latex.
- Glove for professional use.
- In the healthcare environment, this glove is used to make barrier against microorganisms between the healthcare professional, the patient and their environment.
- In an industrial environment, this glove is designed to protect the wearer from tested chemicals and microorganisms (bacteria, molds and viruses).
- The powder-free glove reduces the risk of airborne contamination and limits the risk of irritation.
- For food handling / preparations (HACCP).

#### PRODUCT DESCRIPTION

Product name:	<b>SafeTouch® Vitals</b>
Product type :	Single use, non-sterile
Material :	Nitrile
Colour :	Blue
Shape :	Ambidextrous
Cuff:	Rolled
Surface :	Textured fingertips
Powdered :	No
Length :	240 mm
Palm thickness :	0,07 mm (±0,02mm)
Finger thickness :	0,09 mm (±0,02mm)
Unit weight (M):	3,2 g

Item	Size
1193N-A	T 6 XS
1193N-B	T 7 S
1193N-C	T 8 M
1193N-D	T 9 L
1193N-E	T 10 XL

Marking on packaging:



#### INDUSTRIES

Food / Health / Medical / Dental/ Administration / Hygiene/ Industry

# 1193N NITRILE GLOVE SAFETOUCH VITALS POWDER-FREE

## Class I Medical Device and Category III Personal Protective Equipment

### CERTIFICATION

Manufactured under an ISO 13485 and ISO 9001 certified quality system.  
 In accordance with requirements of Reg (EU) 2017/745, on medical devices.  
 Comply with the standards EN 455-1/2/3/4 for single use medical gloves.  
 In accordance with requirements of Reg (EU) 2016/425, on personal protective equipments.  
 Comply with the EN ISO 21420, EN ISO 374-1/2/4/5, for protective gloves against tested chemicals and microorganisms including viruses.  
 In accordance with requirements of Reg.(EC) 1935/2004 and the French decree of 5th August 2020. This glove is suitable for food contact, with all types of foodstuff.

### PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.  
 The reuse or extended use may lead to infection or cross-contamination.  
 Please refer to your national legislation about disposal of this product once it has been used

### STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored in a cool environment or at room temperature.  
 Shelf life of the product: 5 years.

### LOGISTIC INFORMATION








**Specifications  
outer case**



**Specifications  
inner case**

**Colour code  
inner and outer**

Item	Size mm	Gross weight kg	QTY/pallet	Size mm	QTY	Colour
1193N-A	262 x 252 x 270	3,2	60	245 x 125 x 48	10 x 100	
1193N-B	262 x 252 x 270	3,4	60	245 x 125 x 48	10 x 100	
1193N-C	262 x 252 x 270	3,7	60	245 x 125 x 48	10 x 100	
1193N-D	262 x 252 x 270	4,5	60	245 x 125 x 48	10 x 100	
1193N-E	262 x 252 x 270	5,1	60	245 x 125 x 48	10 x 100	



### EAN Codes

Item	Inner case	Outer case
1193N-A	3 700 381 519 999	3 700 381 519 982
1193N-B	3 700 381 520 018	3 700 381 520 001
1193N-C	3 700 381 520 032	3 700 381 520 025
1193N-D	3 700 381 520 056	3 700 381 520 049
1193N-E	3 700 381 520 070	3 700 381 520 063