

20.5272

MEDICAL BLANKET

White - 150 x 220 cm

FEATURES & BENEFITS

- Used in hospitals, this product is the guarantee of offering a product of excellent quality
- Protects against cold and provides patient comfort
- The raw materials used are soft
- Tear-proof product



PRODUCT DETAILS

Product name:	Medical Blanket
Product type:	Single use, non-sterile
Material:	Polypropylene + Polyester
Confection:	Ultrasonic welding
Colour:	White
Grammage:	50 g/m ² PP + 130 g/m ² Polyester (± 10%)
Size:	150 x 220 cm (± 1 cm)
Unit weight:	595 g (± 10%)
Packaging:	Individual polybag
Origin:	China

INDUSTRIES

Health / Medical

STORAGE CONDITIONS

Normal storage conditions: do not expose to moisture or direct sunlight. Store in a cool environment or at room temperature.

LOGISTICAL INFORMATION

OUTER CASE SPECIFICATIONS					BAG SPECIFICATIONS	
ITEM	DIMENSIONS (mm)	WEIGHT (kg)	Qty / palet	EAN CODE	Packaging	EAN CODE
20.5272	450 x 350 x 420	13	12 boxes (3 layers x 4 boxes)	3 700 381 512 310	20 individuals bags	3 700 381 512 303