

FICHE TECHNIQUE

ADULT BIBS

With pocket

PP + Tissue

66 x 37 cm

Kolmi®



FEATURES & BENEFITS

- Easy to put in place
- One side made of PE : durable and waterproof
- One side made of tissue : absorbent
- One size
- With pocket

TECHNICAL DESCRIPTION

Product name:	Adult bibs with pocket
Product type:	Single use, non sterile
Material :	PE + Tissue
Grammage :	Tissue absorbant : 20 g/m2 Plie PE : 13 g/m2
Unit Weight:	10,4 g (±10%)
Size:	660 x 370 mm (±5mm)
Packaging:	Bag of 100 units
Colour :	White
Origin :	Chine

MARKING ON PACKAGING



* Do not use if package is damaged

LOGISTICAL INFORMATION

OUTER CASE SPECIFICATIONS						
Article	Size (mm)	Weight (kg)	Qty/ palet	Qty/ box	EAN code box	EAN code bag
5605-E	380 x 280 x 550	9,7	30 boxes (5 layers of 6 boxes)	9 bags of 100 units	8 718 007 012 712	8 718 007 912 715