

## 70.498 SAFEWEAR POLYPROPYLENE CAP WITH NECK PROTECTION

### FEATURES & BENEFITS

- The elastic at the back of the head provides excellent support.
- Polypropylene ensures good breathability and thus contributes to the wearer's comfort.
- The neck protector provides good protection against potential contamination.
- The one-size-fits-all design adapts to all body shapes.



### PRODUCT DETAILS

Product name :	Polypropylene cap with neck protection
Product type :	Single use non sterile
Material :	Polypropylene
Confection :	Sewn and elasticated at the back of the head
Dimensions :	32 cm x 30 cm (+/-1cm)
Grammage :	25 g/m <sup>2</sup> (+/-3%)
Unit weight:	3,8 g (+/- 5%)
Color :	White
Qty per case:	5 x 200 pcs
Origin:	China

### CERTIFICATIONS, STANDARDS & TESTS

This product is not a PPE under Regulation 2016/425 governing Personal Protective Equipment.

### STORAGE CONDITIONS

Normal storage and conservation conditions: should not be exposed to moisture and sunlight, should be stored at room temperature.

Product shelf life: 5 years

### LOGISTICS

Item	Qty per case	Case size (cm)	Case gross weight (kg)	Case qty per pallet	Pallet height	EAN Case
70.498	5 x 200 pcs	38 x 28 x 34	4	32 (4 layers of 8 cartons)	1,51 m	3 700 381 522 746