

Technical Data Sheet

PP CLIP CAP – 15g/m² Medical Device Class I

FEATURES & BENEFITS

- 100% Polypropylene spunbonded nonwoven
- Gluing and ultrasound assembly
- LATEX-FREE synthetic elastic under the hem for better comfort
- Soft and very resistant non-woven
- Very breathable for better comfort
- Extensible - suitable for all hair styles



RECOMMENDED USE

- Can be used for nursing staff or for patients in the operating room and in the care departments
- Can also be used in a pharmaceutical, industrial, food industry environments

PRODUCT DETAILS

Product name : Kolmi® – PP Clip Cap 15g/m²
 Product type : Single use, non-sterile
 Material : Polypropylene
 Size : Length : 44,5 cm (± 0,5cm)
 Width : 48 cm (± 1 cm)
 Unit weight : 3,4g (±10%)
 Manufacturing country : FRANCE

| Item | Color |
|--------|-------|
| C12001 | White |
| C12101 | Blue |
| C12201 | Green |

MECHANICAL INFORMATION

Mesurements done by independent laboratories
(measures are available upon request)

| Test | Standard | Laboratory | | Average | Unit | Tolerance |
|------------------|------------|------------|----------------------|---------|------------------|-----------|
| Grammage | ISO 9073-1 | Cetelor | | 15 | g/m ² | ±10% |
| Tensile strength | ISO 9073-3 | Cetelor | Tensile strength SM: | 34,1 | N | ±10% |
| | | | Tensile strength ST: | 15,6 | N | ±10% |
| Elongation | ISO 9073-3 | | Elongation SM: | 44,5 | % | ±10% |
| | | | Elongation ST: | 51,6 | % | ±10% |

PP CLIP CAP - 15g/m²

Medical Device Class I

MICROBIOLOGICAL INFORMATION

Assessment of initial microbial contamination according to standards NF EN ISO 11737: 2018+A1:2021.
Additional microbiological controls: ASR, E.coli, Staphylococci, available upon request.

CERTIFICATIONS & STANDARDS

Meets the requirements of Regulation (EU) 2017/745 on medical devices.
Manufacturing site : ISO 9001 and ISO 13485 certified.


PRECAUTIONS

The product must be used on intact skin only.
The reuse or extended use may lead to infection or cross-contamination.
Please refer to your national legislation on the subject of disposal of this product once it has been used.

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.
Shelf life of the product: 5 years.

MARKING ON PACKAGING



*Do not use if package is damaged

LOGISTICS INFORMATION

| OUTERBOX SPECIFICATIONS | | | | |
|-------------------------|-----------------|-------------|--|-------------------|
| Item | Size (mm) | Weight (kg) | Qty / pallet | Barcode |
| C12001 | 295 x 260 x 190 | 0,9 | 96 cartons (8 layers of 12 cartons) | 3 662 036 000 897 |
| C12101 | | | | 3 662 036 000 910 |
| C12201 | | | | 3 662 036 000 934 |