

Technical Data Sheet

PP CLIP CAP – 15g/m² Medical Device Class I

FEATURES & BENEFITS

- 100% Polypropylene spunbonded nonwoven
- Gluing and ultrasound assembly
- LATEX-FREE synthetic elastic under the hem for better comfort
- Soft and very resistant non-woven
- Very breathable for better comfort
- Extensible - suitable for all hair styles



RECOMMENDED USE

- Can be used for nursing staff or for patients in the operating room and in the care departments
- Can also be used in a pharmaceutical, industrial, food industry environments

PRODUCT DETAILS

Product name : Kolmi® – PP Clip Cap 15g/m²
 Product type : Single use, non-sterile
 Material : Polypropylene
 Size : Length : 44,5 cm (± 0,5cm)
 Width : 48 cm (± 1 cm)
 Unit weight : 3,4g (±10%)
 Manufacturing country : FRANCE

Item	Color
C12501A	White
C12511A	Blue

MECHANICAL INFORMATION

Mesurements done by independent laboratories
(measures are available upon request)

Test	Standard	Laboratory		Average	Unit	Tolerance
Grammage	ISO 9073-1	Cetelor		15	g/m ²	±10%
Tensile strength	ISO 9073-3	Cetelor	Tensile strength SM:	34,1	N	±10%
			Tensile strength ST:	15,6	N	±10%
Elongation	ISO 9073-3		Elongation SM:	44,5	%	±10%
			Elongation ST:	51,6	%	±10%

PP CLIP CAP - 15g/m²

Medical Device Class I

MICROBIOLOGICAL INFORMATION

Assessment of initial microbial contamination according to standards NF EN ISO 11737: 2018+A1:2021.
Additional microbiological controls: ASR, E.coli, Staphylococci, available upon request.

CERTIFICATIONS & STANDARDS

Meets the requirements of Regulation (EU) 2017/745 on medical devices.
Manufacturing site : ISO 9001 and ISO 13485 certified.

PRECAUTIONS

The product must be used on intact skin only.
The reuse or extended use may lead to infection or cross-contamination.
Please refer to your national legislation on the subject of disposal of this product once it has been used.

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.
Shelf life of the product: 5 years.

MARKING ON PACKAGING



*Do not use if package is damaged

LOGISTICS INFORMATION

OUTERBOX SPECIFICATIONS				INNERCASE SPECIFICATIONS			
Item	Size (mm)	Weight (kg)	Qty / pallet	Barcode	Size (mm)	Qty/carton	Barcode
C12501A	800 x 400 x 240	3,4	15 cartons (3 layers of 15 cartons)	3 662 036 000 972	/	12 polybags of 100 unités	3 662 036 003 645
C12511A	800 x 400 x 240	3,4	15 cartons (3 layers of 15 cartons)	3 662 036 000 996	/	12 polybags of 100 unités	3 662 036 003 669