

FICHE TECHNIQUE

WIPES

Spunlace 60 g/m²

28 x 38 cm

Kolmi®



Item	Colour
E09207	☐ White

FEATURES & BENEFITS

- Very absorbent and resistant nonwoven sheet

PRECONISATION:

- Used for wiping body
- Used for surface cleaning

TECHNICAL DESCRIPTION

Product name:	Wipes – Spunlace 60 g/m ²
Product type:	Single use, non sterile
Material :	Spunlaced nonwoven 50% Viscose fibers / 50% Polyester fibers
Grammage :	60 g/m ²
Unit Weight:	6 g (±10%)
Size:	280 x 380 mm (±5mm)
Packaging:	Bag of 70 units
Colour :	White
Origin :	France 

MARKING ON PACKAGING



* Do not use if package is damaged

MECHANICAL INFORMATION

TEST	TEST METHOD		AVERAGE	UNIT	TOLERANCE
Grammage	NWSP 130.1.R0		55	<i>g/m²</i>	<i>± 10%</i>
Thickness	NWSP 120.6.R0		0,57	<i>mm</i>	<i>± 0,12mm</i>
Tensile strength	NWSP 110.4.R0	Dry MD :	85	<i>N/5cm</i>	<i>min 55</i>
		Dry CD :	21	<i>N/5cm</i>	<i>min 14</i>
Elongation		Dry MD :	35	%	<i>max 47</i>
		Dry CD :	125	%	<i>max 155</i>
Absorption	NWSP 010.1.R0		900	%	<i>min 600%</i>

MICROBIOLOGICAL INFORMATION

Assessment of initial microbial contamination in accordance with NF EN ISO 11737: 2018 + A1: 2021.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTICAL INFORMATION

OUTER CASE SPECIFICATIONS					
Article	Size (mm)	Weight (kg)	Qty/ palet	Qty/ box	EAN code bag
E09207	600 x 400 x 219	5,6	32 boxes (4 boxes x 8 layers)	12 bags of 70 units	3 662 036 001 603