

G09010 SAFECARE SPUNLACE WHASH GLOVE 80g/m² 15 x 22cm

FEATURES & BENEFITS

- Absorbent, soft and very resistant washing glove
- Ultrasonic sealing



PRODUCT DESCRIPTION

Product name:	SafeCare Wash glove – Honeycomb spunlace
Product type:	Single use, non-sterile
Material:	Spunlace
Unit Weight:	5g (±10%)
Grammage:	80g/m ²
Dimension:	150 x 220mm (±5mm)
Colour:	White
Packaging:	50 units in bag
Country of origin:	China

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, store in a cool environment or at room temperature.

PRODUCT PERFORMANCES

Test	Standard		Average	Unit	Tolérance
Grammage	NF EN ISO 9073-1		80	g/m ²	± 5%
Thickness	NF EN ISO 9073-2		0,5	mm	± 10%
Tensile Strength	NF EN ISO 9073-3	Dry MD:	169	N/5cm	± 10%
		Dry CD:	119	N/5cm	± 10%
		Wet MD:	136	N/5cm	± 10%
		Wet CD:	124	N/5cm	± 10%
Elongation	NF EN ISO 9073-3	Dry MD:	60	%	± 10%
Absorption	NF EN ISO 9073-6		800	%	± 10%

LOGISTICAL INFORMATION

PACKAGING SPECIFICATIONS

Item	Size (mm)	Weight (kg)	Qty/ Pallet	Outer Case EAN	Qty / Outer Case
G09010	470 x 310 x 360	5,5	20 boxes 4 layers of 5 cartons	3 700 381 511 498	20 Polybags of 50 units