

**H03602A****Blue polypropylene shoecovers without sole**

Medical Device Class I

**FEATURES & BENEFITS**

- Hydrophobic spunbond polypropylene
- Anatomical shape that adapts to all types of shoes and morphologies when used without shoes, particularly in medical settings
- Elastic ankle support

**Recommended uses**

- In industrial environments on dry floors
- In medical environments

**PRODUCT DETAILS**

Product name :	Kolmi®
Product type:	Single-use, non-sterile
Composition:	Spunbond polypropylene Embossed polyethylene film
Spunbond weight :	38 g/m <sup>2</sup> (+/- 10%)
Manufacturing :	Gluing
Size :	190 x 380 mm (Size L + tall height)
Unit weight:	7 g (+/- 10%)
Packaging :	Dispenser case of 500 pcs
Origin :	France

**CERTIFICATIONS & STANDARDS**

Meet requirements of Regulation (EU) 2017/745 on Medical Devices.  
Manufacturing site: ISO 9001 and ISO 13485 certified.

**MICROBIOLOGICAL INFORMATION**

Evaluation of initial microbial contamination according to standards ISO 11737:2018 + A1:2021 upon request.  
Additional microbiological tests: ASR, E. coli, staphylococci, available upon request.

## MECHANICAL RESISTANCE DATAS (Cetelor Laboratory)

Tests on the polypropylene spunbond nonwoven	Standards	Results
Grammage	EN ISO 9073-1	38 g/m <sup>2</sup> (+/-10%)
Tensile strength	EN ISO 9073-3:2023	Machine direction: 83,4 N (+/-10%) Cross direction: 61,08 N (+/-10%)
Elongation at maximum force	EN ISO 9073-3:2023	Machine direction: 53,04 % (+/-10%) Cross direction: 48,08 % (+/-10%)

## NORMAL STORAGE CONDITIONS

Must not be exposed to moisture or sunlight, and should be stored in a cool or room temperature environment.

Product shelf life: 5 years.

## LOGISTICS



Outer case  
spec.

Item	Packaging Type	Qty pcs	Size (mm)	Gross weight (kg)	QTY/pallet
H03602A	Dispenser case	500	400 x 400 x 375	4,1	30

 EAN Codes

Item	Outer case
H03602A	3 662 036 002 501