

1163

LATEX GLOVE SAFETOUCH REJUVENATE POWDER-FREE

Enriched in Lanolin and Vitamin E

Class I Medical Device and Category III Personal Protective Equipment

FEATURES & BENEFITS

- Great elasticity.
- Good dexterity.
- Latex is the material which allows the best contact sensibility.
- Excellent dexterity and grip of objects.
- Enriched in Lanolin and Vitamin E, this glove gives the feeling of a better hydrated skin, softer and smoother.
- Contraindication: Latex is a natural material that contains protein which might cause allergies.

Medicom
SafeTouch®
Rejuvenate

Lanolin &
Vitamin E



Recommendation:

- Glove for professional use.
- In the healthcare environment, this glove is used to make barrier against microorganisms between the professional and the patient and their environment.
- In an industrial environment, this glove is intended to protect the wearer from microorganisms (bacteria, moulds and viruses).
- The powder-free glove reduces the risk of airborne contamination and limits the risk of irritation.

PRODUCT DESCRIPTION

Product name :	SafeTouch Rejuvenate
Product type :	Single use, non-sterile
Material :	Latex
Colour :	Natural
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Textured fingertips
Powdered :	No
Length :	240 mm
Palm thickness :	0,14±0,03 mm
Finger thickness :	0,15±0,03 mm
Integrated lotion:	Lanolin and Vitamin E
Unit weight (M):	5,2g (±0,2g)
Country of origin :	Malaysia

Item	Size	
1163-A	T 6	XS
1163-B	T 7	S
1163-C	T 8	M
1163-D	T 9	L
1163-E	T10	XL

Marking on packaging:



INDUSTRIES

Health / Medical/ Dental

CERTIFICATION

Manufactured under an ISO 13485 and ISO 9001 certified quality system.

In accordance with requirements of Reg (EU) 2017/745, on medical devices.

Comply with the standards EN 455-1/2/3/4 for single use medical gloves.

In accordance with requirements of Reg (EU) 2016/425, on personal protective equipments.

Comply with the EN ISO 21420, EN ISO 374-2/5, ISO 16604 for protective gloves against microorganisms including viruses.

PRECAUTIONS BEFORE USE

The glove must be used on intact skin only.

Single use, do not reuse.

Please refer to your national legislation about disposal of this product once it has been used.

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, store in a cool environment or at room temperature.

Shelf life of the product: 5 years.

LOGISTIC INFORMATION



Specifications outer case

Item	Size mm	Gross weight kg	QTY/ Pallet
1163-A	370 x 255 x 250	6,3	60
1163-B	370 x 255 x 250	6,8	60
1163-C	370 x 255 x 250	7,1	60
1163-D	370 x 255 x 250	8,3	60
1163-E	370 x 255 x 250	8,8	60



Specifications inner case

Size mm	QTY
240 x 120 x 70	10 x 100
240 x 120 x 70	10 x 100
240 x 120 x 70	10 x 100
240 x 120 x 70	10 x 100
240 x 120 x 70	10 x 100

Colour code inner and outer

Colour








EAN Codes

Item	Inner case	Outer case
1163-A	3 700 381 514 970	3 700 381 514 581
1163-B	3 700 381 514 987	3 700 381 514 598
1163-C	3 700 381 514 994	3 700 381 514 604
1163-D	3 700 381 515 007	3 700 381 514 611
1163-E	3 700 381 520 285	3 700 381 520 292